



Weck®  
**Visistat®**  
Skin Stapler

## Feel the Value

The Weck® Visistat® Skin Stapler is designed to meet the veterinary surgeon's need for an easy-to-use instrument with quality performance for external skin closure.

Each disposable Visistat Skin Stapler is loaded with 35 regular or wide medical grade stainless steel staples that can be used in a range of surgical procedures.


The precision instrument features a see-through window designed for better visibility and a patented non-stop feed mechanism designed for easy release and reduced jamming.

The Visistat Skin Stapler is perfectly supplemented by our Weckstat™ Skin Staple Extractor, providing both comfort and simplicity during staple removal.





Staple Extractor ring handles provide improved feel and control while the metal jaws provide strength. Packaged sterile.

# Weck® Visistat® Skin Stapler Features



## Competitive Comparison



## Competitive Comparison


### Staple Height



1 Weck® Visistat®  
35W Skin Staple



2 Dee Veterinary  
Skin Staple 35W



3 Medtronic® Appose®  
ULC Skin Stapler 35W

## Skin Staple Competitive Comparison

PRODUCT NAME	CLOSED STAPLE WIDTH (MM)	CLOSED STAPLE HEIGHT (MM)
1 Weck Visistat 35W Skin Staple	6.5	4.7
2 Dee Veterinary Skin Staple 35W	7.0	4.0
3 Medtronic Appose ULC Skin Stapler 35W	6.5	4.1



Since 1984, RICA Surgical Products has been providing high quality surgical instruments, medical equipment and surgical suite supplies.

www.ricasurgical.com  
1.800.889.3218

Rx only. Not all products may be available in all markets.

RICA is a registered trademark of RICA Surgical Products, Inc.

Medtronic and Appose are registered trademarks of Medtronic.

Competitive Data on File.

Teleflex, the Teleflex logo, Visistat, Weckstat, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 11/2019.

© 2019 Teleflex Incorporated. All rights reserved. MC-005847 Rev 0

Teleflex®

Teleflex • 3015 Carrington Mill Boulevard, Morrisville, NC 27560  
Toll Free: 866 246 6990 | Phone: +1 919 544 8000  
teleflex.com