



PRODUCTS, DEVICES & SYSTEMS

Quality Surgical Suite Products Since 1984



Founded in 1984, RICA Animal Health has been a trusted provider of high-quality surgical instruments and equipment, serving veterinarians, technicians and Animal Health professionals across the United States for over 40 years. RICA Animal Health is a division of RICA Surgical.

Based in Schiller Park, Illinois, RICA Animal Health's foundation is built on a commitment to delivering high quality products at exceptional value to our customers.

As a leading supplier of veterinary surgical instruments, RICA Animal Health offers a wide range of products, including Weck® Visistat® Skin Staplers, Weck® Horizon® Ligation, Paw Protector® Cage Racks, SteriBest® sterilization products, and much more. We also specialize in custom instrument orders, tailored to your exact specifications, and take pride in delivering the right solutions to meet the unique needs of your practice. Contact us today to explore our full line of products.

For pricing on any of the products included in this catalog, please contact us at 1-800-889-3218 or sales@RICAanimalhealth.com.

Thank you for reviewing our product offerings. We appreciate your interest in our surgical suite products. RICA Animal Health will work hard to continue to earn your business by providing you with excellent customer service and a fantastic value for your purchases.

Weck®, Horizon®, Visistat®, Pilling® and KMedic® are registered trademarks of Teleflex Medical, Inc.

TABLE OF CONTENTS

SURGICAL INSTRUMENTATION	4
Instrument Quality Grades	4
RICA Surgical Instrument Guarantee / Warranty	5
Surgical High Temperature Cautery	6
Scissors — Dissecting	7
Scissors — Operating	12
Scissors — Stitch / Suture / Wire Cutting	14
Scissors — Bandage / Universal / Gauze	16
Forceps — Dressing / Tissue	18
Forceps — Hemostatic	23
Needle Holders	30
Knife Handles	35
Towel Clamps / Spay Hooks / Toenail Clippers	37
Tubing Clamps / Sponge Forceps	38
Ophthalmic Instruments	39
Suction	45
Retractors	47
Hooks	54
Orthopedic Instruments	56
OB / GYN Instruments	65
ENT Instruments	67
Additional Instruments	69
Pilling® & KMedic® Instruments	71
Surgical Instrument Packs	76
STERILIZATION PRODUCTS & SYSTEMS	82
Instrument Trays	82
Instrument Enzymatic Cleaning Solution	89
Autoclave Gloves	90
Cleaning Brushes	91
Instrument Cleaning	92
Facial Protection	93
Markers	94
VETERINARY SURGICAL DEVICES & PRODUCTS	95
Weck® Visistat® Skin Stapling System	95
Horizon™ Ligation Clip System	96
BFW™ Veterinary Headlights	98
Firefly Global™ Wireless HD Video Otoscope	99
Mediflex® Dynatrac® Retractor System	100
Tecnimed VisioFocus VET Thermometer	102
Paw Protector® Cage Racks	103
ORDERING INFORMATION / STERIBEST WARRANTY	104
INDEX	105

INSTRUMENT QUALITY GRADES

RICA Animal Health makes it convenient and easy for you to purchase surgical instruments. RICA offers the surgeon a wide variety of surgical instruments with a choice of quality and prices. In addition, RICA fully supports their instruments with a guarantee and unique warranty, as well as a complete instrument repair service. Therefore, you are able to purchase the surgical instruments you need from RICA to meet your needs, whether it be high quality instruments for longevity or economical instruments to stay within your budget.

RICA offers a choice of two grades of surgical instruments varying in quality and cost. These choices allow you to purchase all of the instruments you desire based on your specifications from one manufacturer, simplifying the purchase process and acquisition of your instruments. Make your selection of quality and price from the two grades of surgical instruments RICA offers.

GRADE 1

This grade of surgical instrument is the highest quality instrument available in the market today. These precision German-quality instruments are forged and handcrafted specifically for RICA by master craftsmen. Using the most exacting standards in the industry, these instruments are manufactured using 410-420 surgical grade stainless steel, which is the best series of steel to use for maintaining hardness and longevity. Grade 1 instruments will provide the surgeon with unfaltering quality, excellent durability, consistent reliability and are a great investment.

All instruments in this category include an unconditional lifetime guarantee on materials and workmanship and represent the best value based on long term economics. All inserts for tungsten carbide instruments are warranted for three years.

Instruments of this grade are designated with a catalog number that begins with the prefix number 1 (i.e. 1-xxx). Therefore, if you want to purchase the highest quality grade of surgical instruments, please order the instruments with this number designation.

GRADE 2

The surgical instruments classified as Grade 2 are high quality instruments that are more affordable than our Grade 1 instruments. These instruments are manufactured to meet the rigorous standards and daily requirements of the operating room (hence "operating room" quality). These instruments are manufactured using 410-420 surgical grade stainless steel. The surgical instruments in this grade represent a high quality alternative that will provide more short-term value for the surgeon. An excellent line of instruments to purchase for less critical procedures and when price is important.

These instruments also offer an unconditional lifetime guarantee against defects in materials and workmanship. All inserts for tungsten carbide instruments are warranted for three years.

All instruments of this grade are designated with a catalog number that begins with the prefix number 2 (i.e. 2-xxx). Therefore, if you want to purchase this quality grade of surgical instruments, please order the instruments with this number designation.

RICA SURGICAL INSTRUMENT GUARANTEE / WARRANTY

All RICA Animal Health manufactured surgical instruments are supported by our written lifetime guarantee. RICA's unique 3-year warranty supports all tungsten carbide insert instruments in Grade 1 and Grade 2.

— RICA Animal Health

GUARANTEE

All RICA brand surgical instruments are manufactured to our high standards and have an unconditional lifetime guarantee against defects in materials and workmanship.

Any RICA surgical instrument which proves to be defective in workmanship or materials, will either be repaired or replaced free of charge. This guarantee does not cover normal wear and tear. Misuse, modifications made to any instrument, and/or repair by other than an authorized facility will void this guarantee.

WARRANTY

All RICA brand surgical instruments with tungsten carbide inserts are warranted for three years. RICA will replace these inserts at no charge for the first three years after purchase.

SURGICAL HIGH TEMPERATURE CAUTERY

- Surgical High Temperature Loop Tip Cautery can be used for a wide variety of procedures.
- Cap end protects the delicate loop tip in transport and disposal.
- Cordless
- Packaged sterile



Code	Size
1-6826S	8.5"

SCISSORS — DISSECTING

Code	Description	Size
1-5120		5.5"
2-5120		5.5"
1-5120TC	Carbide Inserts	5.5"
1-5136R	Round Blade	5.5"
1-5124		6.75"
2-5124		6.75"
1-5124TC	Carbide Inserts	6.75"
1-103-5017L	Left Handed	6.75"
1-132-860L	Left Handed, Carbide Inserts	6.75"
1-5128		9.0"
1-5128TC	Carbide Inserts	9.0"
1-5132		10.0"
1-5132TC	Carbide Inserts	10.0"
1-5136		11.0"
1-5136TC	Carbide Inserts	11.0"

MAYO SCISSORS, STRAIGHT



Code	Description	Size
1-5122		5.5"
2-5122		5.5"
1-5122TC	Carbide Inserts	5.5"
1-5118	Round Blade	5.5"
1-5126		6.75"
2-5126		6.75"
1-5126TC	Carbide Inserts	6.75"
1-103-5117L	Left Handed	6.75"
1-132-820L	Left Handed, Carbide Inserts	6.75"
1-5130		9.0"
2-5130		9.0"
1-5130TC	Carbide Inserts	9.0"
1-5134		10.0"
1-5134TC	Carbide Inserts	10.0"
1-5138		11.0"
1-5138TC	Carbide Inserts	11.0"
1-103-5145		18.0"


MAYO SCISSORS, CURVED



SCISSORS — DISSECTING

Code	Description	Size
1-5144	Straight	6.5"
1-5144TC	Straight, Carbide Inserts	6.5"
1-5146	Curved	6.5"
1-5146TC	Curved, Carbide Inserts	6.5"

MAYO NOBLE SCISSORS



Code	Description	Size
1-5160	Straight	7.0"
1-5162	Curved	7.0"

DOYEN FERGUSON SCISSORS



Code	Description	Size
1-271001	Straight	8.0"
1-271000	Curved	8.0"

MARTIN CARTILAGE SCISSORS SERRATED BLADES



Code	Description	Size
1-5220	Straight	8.0"
1-5220TC	Straight, Carbide Inserts	8.0"
1-5222	Curved	8.0"
1-5222TC	Curved, Carbide Inserts	8.0"

SIMS SCISSORS



SCISSORS — DISSECTING

Code	Description	Size
2-181652	Curved	3.75"

**GRADLE SCISSORS, CURVED
RIBBON HANDLE**



Code	Description	Size
1-5177	Straight, Baby Metz	4.5"
1-5179	Straight	5.5"
2-5179	Straight	5.5"
1-5179TC	Straight, Carbide Inserts	5.5"
1-5181	Straight	7.0"
2-5181	Straight	7.0"
1-5181TC	Straight, Carbide Inserts	7.0"
1-5186	Straight	9.0"
2-5186	Straight	9.0"
1-5186TC	Straight, Carbide Inserts	9.0"
1-5190	Straight	11.0"
1-5190TC	Straight, Carbide Inserts	11.0"

METZENBAUM SCISSORS, STRAIGHT



SCISSORS — DISSECTING

Code	Description	Size
1-5178	Curved, <i>Baby Metz</i>	4.5"
1-5180	Curved	5.5"
2-5180	Curved	5.5"
1-5180TC	Curved, Carbide Inserts	5.5"
1-5180ATC	Curved, <i>Delicate, Blunt Tips, Carbide Inserts</i>	5.5"
1-5286	Curved, <i>Delicate</i>	5.5"
1-104-3514L	Curved, <i>Delicate, Left Handed</i>	5.75"
1-132-711L	Curved, <i>Delicate, Left Handed, Carbide Inserts</i>	6.75"
1-5182	Curved	7.0"
2-5182	Curved	7.0"
1-5182TC	Curved, Carbide Inserts	7.0"
1-132-750L	Curved, <i>Left Handed, Carbide Inserts</i>	7.0"
1-132-751L	Curved, <i>Delicate, Left Handed, Carbide Inserts</i>	7.0"
1-5184	Curved	8.0"
1-5184TC	Curved, Carbide Inserts	8.0"
1-5188	Curved	9.0"
2-5188	Curved	9.0"
1-5188TC	Curved, Carbide Inserts	9.0"
1-5192	Curved	11.0"
2-5192	Curved	11.0"
1-5192TC	Curved, Carbide Inserts	11.0"

METZENBAUM SCISSORS, CURVED



Code	Description	Size
1-5312	Straight	4.0"
1-5314	Curved	4.0"


STRABISMUS SCISSORS



SCISSORS — DISSECTING

Code	Description	Size
1-5300	Straight	4.0"
2-5300	Straight	4.0"
1-5304	Straight	4.5"
1-5304	Straight	4.5"
1-5304TC	Straight, Carbide Inserts	4.5"
1-5304LTC	Straight, Left Handed, Carbide Inserts	4.5"

IRIS SCISSORS, STRAIGHT



Code	Description	Size
1-5301	Curved, Delicate	3.5"
1-5302	Curved	4.0"
2-5302	Curved	4.0"
1-5306	Curved	4.5"
2-5306	Curved	4.5"
1-5306TC	Curved, Carbide Inserts	4.5"
1-5306LTC	Curved, Left Handed, Carbide Inserts	4.5"

IRIS SCISSORS, CURVED



Code	Description	Size
1-181462	Straight	4.2"
1-181466	Curved	4.2"
2-181466	Curved	4.2"


STEVENS TENOTOMY SCISSORS



SCISSORS — OPERATING

Code	Description	Size
1-512	Straight	4.5"
2-512	Straight	4.5"
1-514	Straight	5.0"
1-516	Straight	5.5"
2-516	Straight	5.5"
1-516TC	Straight, <i>Carbide Inserts</i>	5.5"
1-518	Straight	6.5"

OPERATING SCISSORS, STRAIGHT Sharp - Blunt Points



Code	Description	Size
1-52	Straight	4.5"
1-54	Straight	5.0"
2-54	Straight	5.0"
1-56	Straight	5.5"
2-56	Straight	5.5"
1-56TC	Straight, <i>Carbide Inserts</i>	5.5"

OPERATING SCISSORS, STRAIGHT Sharp - Sharp Points



Code	Description	Size
1-522	Straight	4.5"
1-524	Straight	5.0"
2-524	Straight	5.0"
1-526	Straight	5.5"
2-526	Straight	5.5"
1-526TC	Straight, <i>Carbide Inserts</i>	5.5"

OPERATING SCISSORS, STRAIGHT Blunt - Blunt Points



SCISSORS — OPERATING

Code	Description	Size
1-542	Curved	4.5"
1-544	Curved	5.0"
2-544	Curved	5.0"
1-546	Curved	5.5"
2-546	Curved	5.5"
1-546TC	Curved, Carbide Inserts	5.5"
1-549	Curved	7.0"

OPERATING SCISSORS, CURVED Sharp - Blunt Points



Code	Description	Size
1-532	Curved	4.5"
2-532	Curved	4.5"
1-534	Curved	5.0"
1-536	Curved	5.5"
1-536TC	Curved, Carbide Inserts	5.5"
1-539	Curved	7.25"

OPERATING SCISSORS, CURVED Sharp - Sharp Points



Code	Description	Size
1-552	Curved	4.5"
2-552	Curved	4.5"
1-554	Curved	5.0"
1-556	Curved	5.5"
1-556TC	Curved, Carbide Inserts	5.5"

OPERATING SCISSORS, CURVED Blunt - Blunt Points



SCISSORS — STITCH / SUTURE / WIRE CUTTING

Code	Size
1-9105	4.5"
1-9104	5.5"
2-9104	5.5"

LITTAUER STITCH SCISSORS



Code	Size
1-9100	3.5"
2-9100	3.5"
1-9102	5.0"

SPENCER STITCH SCISSORS



Code	Size
1-9103	4.75"
2-9103	4.75"

NORTHBENT STITCH SCISSORS



Code	Size
1-9101	3.5"
2-9101	3.5"

SHORTBENT STITCH SCISSORS



SCISSORS — STITCH / SUTURE / WIRE CUTTING

Code	Size
1-996	6.25"

HEATH SCISSORS



Code	Size
1-9124	4.75"
2-9124	4.75"

ROGER WIRE SCISSORS



Code	Description	Size
1-9120	Straight	4.75"
1-9122	Curved	4.75"

CROWN & COLLAR SCISSORS ONE SERRATED BLADE



SCISSORS — BANDAGE / UNIVERSAL / GAUZE

Code	Description	Size
1-5510		3.5"
1-5512		4.5"
2-5512		4.5"
1-5514		5.5"
2-5514		5.5"
1-5516		7.5"
2-5516		7.5"
1-5516LH	Left Handed	7.5"

LISTER BANDAGE SCISSORS



Code	Size
1-5550	8.0"
2-5550	8.0"

LISTER BANDAGE SCISSORS ONE LARGE FINGER RING



Code	Size
1-5998	6.5"
1-51000	7.5"

UNIVERSAL (EMT) SHEARS



SCISSORS — BANDAGE / UNIVERSAL / GAUZE

Code	Description	Size
1-5560	Straight	5.5"
2-5560	Straight	5.5"

KNOWLES BANDAGE SCISSORS



Code	Description	Size
1-5560	Angled	5.5"

KNOWLES BANDAGE SCISSORS ONE SERRATED BLADE



Code	Size
1-5578	7.0"

SCHMIEDEN-TAYLOR FELINE BANDAGE SCISSORS



Code	Size
1-5564	8.0"
2-5564	8.0"

U.S.A.. GAUZE SHEARS



FORCEPS - DRESSING / TISSUE

Code	Size
1-64	4.5"
1-66	5.0"
2-66	5.0"
1-610	6.0"
2-610	6.0"
1-614	8.0"
2-614	8.0"
1-616	10.0"
2-616	10.0"
1-617	12.0"

DRESSING FORCEPS - SERRATED TIPS



Code	Size
1-640	4.5"
1-642	5.0"
2-642	5.0"
1-644	5.5"
2-644	5.5"
1-646	6.0"
2-646	6.0"
1-648	8.0"
2-648	8.0"
1-650	10.0"
2-650	10.0"
1-106-0545	18.0"

TISSUE FORCEPS - 1 X 2 TEETH



Code	Size
1-660	4.5"
1-662	5.0"
2-662	5.0"
1-664	5.5"
2-664	5.5"
1-666	6.0"

TISSUE FORCEPS - 2 X 3 TEETH



FORCEPS – DRESSING / TISSUE

Code	Description	Size
1-6124		4.75"
2-6124		4.75"
1-6124XL	<i>Fenestrated Handle</i>	4.75"
1-6126	<i>Angled</i>	4.75"

BROWN ADSON TISSUE FORCEPS - 7 X 7 TEETH



Code	Size
1-6125	4.75", 1.5 mm

DEBAKEY ADSON ATRAUMATIC TISSUE FORCEPS



Code	Size
1-24573	6.0", 1.0 mm
1-24570	6.0", 1.5 mm
1-24571	6.0", 2.0 mm
1-24550	6.0", 2.5 mm
1-24572	8.0", 1.5 mm

DEBAKEY ATRAUMATIC TISSUE FORCEPS



FORCEPS- DRESSING / TISSUE

Code	Description	Size
1-6118		4.75"
2-6118		4.75"
1-6118D	Delicate	4.75"
1-172510	Micro	4.75"

ADSON DRESSING FORCEPS - SERRATED TIPS



Code	Description	Size
1-6120		4.75"
2-6120		4.75"
1-6120L	Long Pattern	4.75"
1-6120XL	Fenestrated Handle	4.75"
1-6123XL	With Tying Platform	4.75"
1-172500	Micro	4.75"

ADSON TISSUE FORCEPS - 1 X 2 TEETH



FORCEPS — DRESSING / TISSUE

Code	Size
#2 - 17302	4.75"
#3 - 17303	4.75"
#4 - 17304	4.75"
#5 - 17305	4.75"
#6 - 17306	4.75"
#7 - 17307	4.75"

JEWELER FORCEPS



Code	Size
1-6154	5.0"

BEASLEY BABCOCK FORCEPS - 6 X 6 TEETH



Code	Size
1-6104	5.0"

SEMKEN DRESSING FORCEPS



Code	Description	Size
1-6106	Straight	5.0"
1-6107	Curved	5.0"

SEMKEN TISSUE FORCEPS - 1 X 2 TEETH



FORCEPS — DRESSING / TISSUE

Code	Size
1-19371	6.25"

JANSEN DRESSING FORCEPS



Code	Size
1-6142	6.0"
2-6142	6.0"
1-6144	8.0"
1-6145	10.0"

RUSSIAN TISSUE FORCEPS



Code	Description	Size
1-6312	Straight	4.25"
1-6314	Curved	4.25"

WALTER SPLINTER FORCEPS



FORCEPS – HEMOSTATIC

Code	Description	Size
1-724	Straight	3.5"
2-724	Straight	3.5"
1-724D	Straight, Delicate	3.5"
1-726	Curved	3.5"
2-726	Curved	3.5"
1-726D	Curved, Delicate	3.5"

HARTMAN MOSQUITO FORCEPS



Code	Description	Size
1-72	Straight	5.0"
2-72	Straight	5.0"
1-78	Straight, Micro	5.0"
1-7250	Straight	8.0"
1-74	Curved	5.0"
2-74	Curved	5.0"
1-710	Curved, Delicate	5.0"
1-74M	Curved, Micro	5.0"
1-7252	Curved	8.25"
1-7254	Right Angle	7.75"

HALSTED MOSQUITO FORCEPS



Code	Description	Size
1-714	Straight	5.0"
2-714	Straight	5.0"
1-716	Curved	5.0"
2-716	Curved	5.0"

HALSTED MOSQUITO FORCEPS - 1 X 2 TEETH



FORCEPS – HEMOSTATIC

Code	Description	Size
1-736	Straight	5.5"
2-736	Straight	5.5"
1-738	Curved	3.5"
2-738	Curved	3.5"

KELLY FORCEPS



Code	Description	Size
1-742	Straight	5.5"
2-742	Straight	5.5"
1-746	Straight	6.5"
2-746	Straight	6.5"
1-744	Curved	5.5"
2-744	Curved	5.5"
1-748	Curved	6.5"
2-748	Curved	6.5"

CRILE FORCEPS



Code	Description	Size
1-7170	Straight	6.25"
2-7170	Straight	6.25"
1-7180	Straight	8.0"
2-7180	Straight	8.0"
1-1163710	Straight	9.0"
1-1163720	Straight	10.0"
1-7172	Curved	6.25"
2-7172	Curved	6.25"
1-7182	Curved	8.0"
2-7182	Curved	8.0"
1-1163725	Curved	10.0"
1-1163730	Curved	12.0"

ROCHESTER CARMALT FORCEPS



FORCEPS — HEMOSTATIC

Code	Description	Size
1-7260	Straight	6.5"
1-7264	Straight	7.5"
1-7262	Curved	6.5"
1-7266	Curved	7.5"

FERGUSON ANGIOTRIBE FORCEPS



Code	Description	Size
1-770	Straight, "Baby"	5.5"
1-7150	Straight	6.25"
2-7150	Straight	6.25"
1-7152	Straight	7.25"
2-7152	Straight	7.25"
1-7154	Straight	8.0"
2-7154	Straight	8.0"
1-7156	Straight	9.0"
1-7157	Straight	10.0"
2-7157	Straight	10.0"
1-7160	Curved	6.25"
2-7160	Curved	6.25"
1-7162	Curved	7.25"
2-7162	Curved	7.25"
1-7164	Curved	8.0"
2-7164	Curved	8.0"
1-7166	Curved	9.0"
1-7167	Curved	10.0"
2-7167	Curved	10.0"

ROCHESTER OSCHNER FORCEPS - 1 X 2 TEETH



FORCEPS — HEMOSTATIC

Code	Description	Size
1-792	Straight	6.25"
1-794	Curved	6.25"
2-794	Curved	6.25"

RANKIN FORCEPS



Code	Description	Size
1-7118	Straight	5.5"
1-7120	Straight	6.25"
1-7122	Straight	7.25"
1-7124	Straight	8.0"
1-7126	Straight	9.0"
1-7128	Straight	10.25"
1-7130	Straight	12.0"
1-7132	Straight	14.0"
1-7136	Curved	5.5"
2-7136	Curved	5.5"
1-7138	Curved	6.25"
2-7138	Curved	6.25"
1-7140	Curved	7.25"
1-7142	Curved	8.0"
2-7142	Curved	8.0"
1-7144	Curved	9.0"
1-7146	Curved	11.0"
1-7148	Curved	12.0"

ROCHESTER PEAN FORCEPS



FORCEPS — HEMOSTATIC

Code	Description	Size
1-16109	Straight, "Baby"	7.0"
1-16160	Straight	9.0"
2-16160	Straight	9.0"
1-16110	Curved, "Baby"	7.0"
1-16162	Curved	9.0"

DOYEN FORCEPS



Code	Size
1-243261	6.5"

COOLEY-SATINSKY FORCEPS



Code	Description	Size
1-163	4 x 5 Teeth, Delicate	5.0"
1-164	4 x 5 Teeth, Delicate	5.5"
1-166	4 x 5 Teeth	6.0"
2-166	4 x 5 Teeth	6.0"
1-168	5 x 6 Teeth	6.0"
1-1610	5 x 6 Teeth	7.5"
1-1611	5 x 6 Teeth	9.0"
1-1612	5 x 6 Teeth	10.0"
1-1628	5 x 6 Teeth, Angled Jaw	12.0"

ALLIS TISSUE FORCEPS



FORCEPS — HEMOSTATIC

Code	Size
1-1644	6.25", 9.0 mm Jaw
1-1646	8.25", 10 mm Jaw
1-1648	9.25", 15 mm Jaw

BABCOCK INTESTINAL FORCEPS



Code	Description	Size
1-7210	Fully Curved, "Baby", Delicate	5.5"
1-7200	Fully Curved	6.25"
2-7200	Fully Curved	6.25"
1-7202	Fully Curved	7.25"
1-7206	Fully Curved, "Baby"	9.0"
1-7212	Right Angle Jaws, "Baby"	5.5"
1-25810	Right Angle Jaws	7.25"
1-25812	Right Angle Jaws	9.0"

MIXTER FORCEPS



Code	Description	Size
1-24804	90° Angle	4.75"

DEBAKEY RING HANDLED BULLDOG CLAMP



FORCEPS — HEMOSTATIC

Code	Dimensions
1-497230	7 mm x 10 mm Wide 9 mm x 10 mm Wide
1-497231	11 mm x 18 mm Wide 12 mm x 18 mm Wide
1-497232	13 mm x 25 mm Wide 14 mm x 25 mm Wide

FERGUSON GALLSTONE SCOOPS



Code	Size
1-1699	9.0"

DENNIS ANASTOMOSIS CLAMP



Code	Size
1-24703	2.0", 3 cm Jaw

GLOVER BULLDOG CLAMP



NEEDLE HOLDERS

Code	Description	Size
1-842		5.5"
1-842TC	Carbide Inserts	5.5"
1-844		6.0"
2-844		6.0"
1-844TC	Carbide Inserts	6.0"
1-844L	Left Handed	6.0"
1-844LTC	Left Handed, Carbide Inserts	6.0"
1-846		7.0"
2-846		7.0"
1-846TC	Carbide Inserts	7.0"
1-848		8.0"
2-848		8.0"
1-848TC	Carbide Inserts	8.0"
1-8130TC	Masson Needle Holder, Carbide Inserts	10.5"

MAYO-HEGAR NEEDLE HOLDER



Code	Description	Size
1-814	<i>Serrated Jaws</i>	4.75"
1-814TC	<i>Serrated Jaws, Carbide Inserts</i>	4.75"
2-814TC	<i>Serrated Jaws, Carbide Inserts</i>	4.75"
1-814ATC	<i>Smooth Jaws, Carbide Inserts</i>	4.75"
1-814LTC	<i>Serrated Jaws, Left Handed, Carbide Inserts</i>	4.75"
1-815		5.5"
2-815		5.5"
1-815TC	<i>Carbide Inserts</i>	5.5"
2-815TC	<i>Carbide Inserts</i>	5.5"
1-815LTC	<i>Left Handed, Carbide Inserts</i>	5.5"
1-816		6.5"
2-816		6.5"
1-816TC	<i>Carbide Inserts</i>	6.5"
2-816TC	<i>Carbide Inserts</i>	6.5"
1-816LTC	<i>Left Handed, Carbide Inserts</i>	6.5"
1-817		7.5"
2-817		7.5"
1-817TC	<i>Carbide Inserts</i>	7.5"
2-817TC	<i>Carbide Inserts</i>	7.5"

OLSEN-HEGAR NEEDLE HOLDER



NEEDLE HOLDERS

Code	Description	Size
1-828	Spring Handle, Serrated Jaws	5.5"
1-828TC	Serrated Jaws, Carbide Inserts	5.5"
1-830	Spring Handle, Serrated Jaws	6.25"
1-830TC	Serrated Jaws, Carbide Inserts	6.25"

MATHIEU NEEDLE HOLDER



Code	Description	Size
1-833TC	Carbide Inserts	5.5"
1-834TC	Carbide Inserts	6.25"
1-835TC	Carbide Inserts	8.0"

MATHIEU OLSEN-HEGAR NEEDLE HOLDER

Code	Description	Size
1-850		6.0"
1-850TC	Carbide Inserts	6.0"
1-8750TCLH	Left Handed, Carbide Inserts	6.0"
1-851		7.0"
1-851TC	Carbide Inserts	7.0"
1-852TC	Carbide Inserts	8.0"
1-132-313TC	Carbide Inserts	18.0"

CRILE WOOD NEEDLE HOLDER



Code	Size
1-822	5.25"
2-822	5.25"

BROWN NEEDLE HOLDER



NEEDLE HOLDERS

Code	Description	Size
1-896		7.2"
1-896TC	Carbide Inserts	7.2"

SAROT NEEDLE HOLDER



Code	Description	Size
1-87ATC	Delicate, Carbide Inserts	4.5"
1-86		5.0"
1-86TC	Carbide Inserts	5.0"

WEBSTER NEEDLE HOLDER, SMOOTH JAW



Code	Description	Size
1-84		4.75"
2-84		4.75"
1-84TC	Carbide Inserts	4.75"

DERF NEEDLE HOLDER



Code	Description	Size
1-88		5.0"
1-88TC	Carbide Inserts	5.0"


HALSEY NEEDLE HOLDER



NEEDLE HOLDERS


Code	Description	Size
2-8688TC	Carbide Inserts	5.25"

SIDE CUTTING NEEDLE HOLDER



Code	Description	Size
1-870		7.0"
1-870TC	Carbide Inserts	7.0"

ADSON NEEDLE HOLDER



Code	Description	Size
1-181820	Straight, Without Lock	5.5"
1-181828	Straight, With Lock	5.5"
2-181828	Straight, With Lock	5.5"
1-181828TC	Straight, With Lock, Carbide Inserts	5.5"
1-181824	Curved, Without Lock	5.5"
1-181832	Curved, With Lock	5.5"
1-181832TC	Curved, With Lock, Carbide Inserts	5.5"

CASTROVIEJO NEEDLE HOLDER



NEEDLE HOLDERS

Code	Description	Size
1-181834	Straight, Without Lock	5.0"
1-181836	Straight, With Lock	5.0"
1-181838	Curved, Without Lock	5.0"
2-181838	Curved, Without Lock	5.0"
1-181839	Curved, With Lock	5.0"

BARRAQUER NEEDLE HOLDER



Code	Description	Size
1-8101TC	Carbide Inserts	5.0"
1-8103TC	Carbide Inserts	6.0"
1-8104TC	Carbide Inserts	7.0"
1-8105TC	Carbide Inserts	8.0"
1-8106TC	Carbide Inserts	9.0"
1-8107TC	Carbide Inserts	10.25"
1-866TC	1.0 mm Jaws, Carbide Inserts	6.0"
1-869TC	Delicate, Carbide Inserts	9.0"

RYDER NEEDLE HOLDER



KNIFE HANDLES

Code	Size
1-47	5.0"
2-47	5.0"

KNIFE HANDLE #3



Code	Size
1-48	5.375"
2-48	5.375"

KNIFE HANDLE #4



Code	Description	Size
1-49	Fits Blades 10-15	6.5"
2-49	Fits Blades 10-15	6.5"

KNIFE HANDLE #7



Code	Description	Size
1-416	Fits Blades 10-15	4.5"
2-416	Fits Blades 10-15	4.5"

KNIFE HANDLE #9



KNIFE HANDLES

Code	Description	Size
1-410	Fits Blades 10-15	8.375"
2-410	Fits Blades 10-15	8.375"

KNIFE HANDLE #3L



Code	Description	Size
1-411	Fits Blades 20-25	8.375"
2-411	Fits Blades 20-25	8.375"

KNIFE HANDLE #4L



Code	Description	Size
1-417	Fits Blades 10-15	5.75"

SIEGEL CORWIN #3 KNIFE HANDLE



TOWEL CLAMPS / SPAY HOOKS / TOENAIL CLIPPERS

Code	Size
1-7504	3.5"
2-7504	3.5"
1-7506	3.5"
2-7506	3.5"

BACKHAUS TOWEL CLAMP



Code	Size
1-7501	2.25"
1-7500	3.5"
2-7500	3.5"

JONES TOWEL CLAMP



Code	Size
1-7514	3.5"
1-7515	5.25"
2-7515	5.25"

LORNA TOWEL CLAMP - NON-PERFORATING



Code	Description	Size
1-28131		8.0"
2-28131		8.0"
1-28131N	Narrow Blade	8.0"
1-28132WB	Wide Blade, Flat Pattern	8.0"
1-28133	Covault Spay Hook	8.0"

SPAY HOOK



Code	Size
1-1718	4.5"

WHITE TOENAIL CLIPPER



TUBING CLAMPS / SPONGE FORCEPS

Code	Size
1-97590	7.0"

PRESBYTERIAN TUBING CLAMP



Code	Description	Size
1-7604	Straight	7.0"
1-7600	Straight	9.5"
2-7600	Straight	9.5"
1-7640	Straight	14.0"
1-7614	Curved	7.0"
1-7610	Curved	9.5"
2-7610	Curved	9.5"

FOERSTER SPONGE FORCEPS - SERRATED JAWS



OPHTHALMIC INSTRUMENTS

Code
1-18650

Size
3.25"

CASTROVIEJO CALIPER GRADUATED IN 0-20 mm



Code
1-181202
2-181202

Inside Ring Size
12 mm X 23 mm
12 mm X 23 mm

Size
3.5"
3.5"

DESMARRES CHALAZION FORCEPS MEDIUM SIZE



Code
1-18833

Size
4.0"

CASTROVIEJO COLIBRI CORNEAL FORCEPS WITH TYING PLATFORM



Code
1-18866

Width
0.5 mm

Size
3.4"

BISHOP-HARMON FORCEPS SERRATED JAWS



OPHTHALMIC INSTRUMENTS

Code	Width	Size
1-18867	0.3 mm	3.4"
2-18867	0.3 mm	3.4"
1-18869	0.5 mm	3.4"
1-18868	0.7 mm	3.4"

BISHOP-HARMON FORCEPS 1 X 2 TEETH



Code	Description	Size
1-18785	Straight, Delicate	4.0"
1-18786	Straight	4.0"
1-18788	Half Curve	4.0"

IRIS FORCEPS 1 X 2 TEETH



Code	Description	Size
1-18780	Straight	4.0"
1-18782	Half Curve	4.0"
1-18784	Full Curve	4.0"

IRIS DRESSING FORCEPS



Code	Width	Description	Size
1-357647	0.5 mm	Round Handle	7.0"
1-357645	1.5 mm	Round Handle	7.0"
1-357681TC	1.0 mm	Round Handle, Carbide Inserts	8.25"

MICRO RING TIP FORCEPS



OPHTHALMIC INSTRUMENTS


Code	Width	Size
1-18827	0.4 mm	3.0"

**TROUTMAN FORCEPS 1 X 2 TEETH
WITH TYING PLATFORM**



Code	Width	Size
1-18950	0.9 mm	3.0"
1-18951	0.3 mm	4.0"

**CASTROVIEJO FORCEPS 1 X 2 TEETH
WITH TYING PLATFORM**



Code	Description
1-1812	15 mm X 5 mm Blades
1-1814	13 mm X 5 mm Blades

CASTROVIEJO EYE SPECULUM



Code	Description	Size
1-1835	Large Blades (15 mm)	1.75"
2-1835	Large Blades (15 mm)	1.75"
1-1836	Small Blades (10 mm)	1.75"
2-1836	Small Blades (10 mm)	1.75"

BARRAQUER WIRE SPECULUM



OPHTHALMIC INSTRUMENTS

Code	Description	Size
1-181620	Straight	3.5"
1-181622	Curved	3.5"
1-181624	Angled on Flat	3.5"

VANNAS SCISSORS




Code	Description	Size
1-394441	Straight, Pencil Grip	6.75"
1-39442	Curved, Pencil Grip	7.0"

TEWES BARRAQUER MICRO DISSECTING SCISSORS



Code	Size
1-181576	3.75"

CASTROVIEJO SCISSORS CURVED



Code	Description	Size
1-181510	Straight	4.5"
1-181516	Curved	4.5"

NOYES IRIS SCISSORS SHARP / SHARP



OPHTHALMIC INSTRUMENTS

Code	Description	Size
1-18450	Small	5.5"
1-18452	Medium	5.5"
1-18454	Large	5.5"

VON GRAEFE STRABISMUS HOOK



Code	Description	Size
1-181820	Straight, Without Lock	5.5"
1-181828	Straight, With Lock	5.5"
2-181828	Straight, With Lock	5.5"
1-181828TC	Straight, With Lock, Carbide Inserts	5.5"
1-181824	Curved, Without Lock	5.5"
1-181832	Curved, With Lock	5.5"
1-181832TC	Curved, With Lock, Carbide Inserts	5.5"

CASTROVIEJO NEEDLE HOLDER



Code	Description	Size
1-181834	Straight, Without Lock	5.5"
1-181836	Straight, With Lock	5.5"
1-181838	Curved, Without Lock	5.5"
2-181838	Curved, Without Lock	5.5"
1-181839	Curved, With Lock	5.5"

BARRAQUER NEEDLE HOLDER



Code	Description	Size
1-18110	Size 1	11 mm
1-18112	Size 2	13 mm
1-18114	Size 3	15 mm
1-18116	Size 4	17 mm

DESMARRES LID RETRACTOR



OPHTHALMIC INSTRUMENTS

Code	Description	Size
1-5150027	Slight Curve	27 Gauge

AIR INJECTION CANNULA



Code	Description	Size
1-1872	Straight	23 Gauge
1-1874	Curved	23 Gauge

LACRIMAL CANNULA



Code	Description	Size
1-18499	Size 0	5.0"
1-18500	Size 1	5.0"
1-18502	Size 2	5.0"
1-18504	Size 3	5.0"
1-18506	Size 4	5.0"

MEYERHOEFER CHALAZION CURETTE



SUCTION


Code	Size
1-26782	8.0", 15 French

ADSON SUCTION TUBE




Code	Size
1-19580	3 French
1-19582	5 French
1-19584	7 French

BARON SUCTION TUBE




Code	Size
1-2104	11.75"

YANKAUER SUCTION TUBE



Code	Size
1-2110	10"

ANDREWS-PYNCHON SUCTION TUBE



SUCTION


Code	Size
1-19570	6 French
1-19572	7 French
1-19573	8 French
1-19574	9 French
1-19575	10 French
1-19576	11 French
1-19577	12 French

FRAZIER-FERGUSON SUCTION TUBE




Code	Description	Size
1-10310	Curved	23 French
2-10310	Curved	23 French
1-10312	Straight	30 French

POOLE SUCTION TUBE



FRENCH CATHETER SCALE



RETRACTORS

Code	Teeth	Size
1-11600	2 x 3 Teeth	4.5"
1-11608	3 x 4 Teeth	5.5"
1-11614	3 x 4 Teeth	6.5"
1-11617	3 x 4 Teeth	8.0"

WEITLANER RETRACTOR SHARP BLADES




Code	Teeth	Size
1-11602	2 x 3 Teeth	4.5"
1-11610	3 x 4 Teeth	5.5"
1-11616	3 x 4 Teeth	6.5"
1-11618	3 x 4 Teeth	8.0"

WEITLANER RETRACTOR BLUNT BLADES




Code	Size
1-11364	3.5"
2-11364	3.5"
1-11362	4.5"
1-11365	5.5"
2-11365	5.5"
1-11366	6.75"
2-11366	6.75"
1-11367	7.5"

GELPI RETRACTOR



Code	Size
1-11368	7.0"

GELPI RETRACTOR WITH BALL STOPS



RETRACTORS

Code	Description	Size
1-1112	Spread of 2.25", Prongs 1/8" Deep	2.75"
2-1112	Spread of 2.25", Prongs 1/8" Deep	2.75"
1-1113	Spread of 3.0", Prongs 1/2" Deep	4.0"

Code	Teeth	Size
1-116	4 x 4 Teeth	4.0"


Code	Teeth	Size
1-112	3 x 3 Teeth	4.25"

Code	Description	Size
1-1114	Sharp Prongs	4.0"
1-1116	Blunt Prongs	4.0"


ALM RETRACTOR



ALLPORT RETRACTOR SHARP PRONGS



JANSEN RETRACTOR SHARP BLADES 3.0" SPREAD



HEISS SELF RETAINING RETRACTOR



RETRACTORS

Code	Description	Size
1-11430	Pediatric, "Baby"	3.5" Spread
1-11420		7.0" Spread
2-11420		7.0" Spread
1-11421		10.0" Spread

BALFOUR ABDOMINAL RETRACTOR



Code	Description	Size
1-1174	Sharp	6.25"
1-1176	Blunt	6.25"
2-1176	Blunt	6.25"

SENN RETRACTOR



Code	Size
1-11114	5.0"
1-11116	7.0"

PARKER RETRACTOR (SET OF 2)



Code	Size
1-11122	8.5"
2-11122	8.5"

U.S. ARMY RETRACTOR (SET OF 2)



RETRACTORS

Code	Blade Size	Size
1-11180	1.0" x 1.75"	10.5"
1-11182	1.125" x 2.875"	10.5"
1-11184	1.125" x 3.25"	10.5"

HEANEY RETRACTOR




Code	Size
1-1173	6.0"
2-1173	6.0"

RAGNELL RETRACTOR



Code	Description	Size
1-11502	1 Prong	8.5"
1-11530	2 Prongs	8.5"
1-11532	3 Prongs	8.5"
1-11534	4 Prongs	8.5"
1-11536	6 Prongs	8.5"

VOLKMAN RETRACTOR SHARP



Code	Description	Size
1-11503	1 Prong	8.5"
1-11538	2 Prongs	8.5"
1-11540	3 Prongs	8.5"
1-11542	4 Prongs	8.5"
1-11544	6 Prongs	8.5"

VOLKMAN RETRACTOR BLUNT



RETRACTORS

Code	Blade Size	Size
1-27213	6 mm	6.25"
1-27214	8 mm	6.25"

MINI-HOHMANN RETRACTOR



Code	Blade Size	Size
1-27215	10 mm	9.25"
1-27216	18 mm	9.25"
1-27217	25 mm	9.25"
1-27218	43 mm	9.25"

HOHMANN RETRACTOR



Code	Blade Size	Size
1-11250	1.0" x 0.75"	9.5"
1-11252	1.25" x 1.0"	9.5"
1-11254	1.5" x 1.0"	9.5"
1-11256	2.0" x 0.75"	9.5"

RICHARDSON RETRACTOR



Code
1-11260

RICHARDSON-EASTMAN RETRACTOR (SET OF 2)



RETRACTORS

Code	Size
1-11112	4.75"
2-11112	4.75"

FARABEUF RETRACTOR (SET OF 2)



Code	Teeth	Size
1-11700	4 x 4 Teeth	7.5"

ADSON RETRACTOR SHARP



Code	Blade Size	Size
1-11200	19 mm	7.0"
1-11202	22 mm	8.5"
1-11204	25 mm	9.0"
1-11206	25 mm	12.0"
1-11208	25 mm	13.0"
1-11210	25 mm	14.0"
1-11212	38 mm	12.0"
1-11216	75 mm	12.0"

DEAVER RETRACTOR



RETRACTORS

Code **Size**
1-11102 9.0"

VEIN (CUSHING) RETRACTOR



Code **Size**
1-18102 5.25"

KNAPP RETRACTOR, 4 BLUNT PRONGS



Code **Size**
1-18100 5.25"

ROLLET RETRACTOR, 4 SHARP PRONGS



HOOKS

Code	Size
1-11138	7.0"

MALLEABLE RIBBON RETRACTOR



Code	Size
1-1106535C	5.0"

WALLACE STIFLE RETRACTOR



Code	Description	Size
1-21153	One Sharp Prong	6.25"
1-21154	Two Sharp Prongs, 2 mm	6.25"
1-21156	Two Sharp Prongs, 5 mm	6.25"
1-21158	Two Sharp Prongs, 7 mm	6.25"
1-21160	Two Sharp Prongs, 10 mm	6.25"

JOSEPH HOOK



HOOKS

Code	Description	Size
1-2185	Small	5.0"
1-2186	Medium	5.0"
1-2187	Large	5.0"

KLEINERT-KUTZ HOOK



Code	Description	Size
1-18416	Sharp Double Prongs	5.0"

GUTHRIE FIXATION HOOK



Code	Description	Size
1-261084	Sharp Point	5.0"

FRAZIER DURA HOOK



ORTHOPEDIC INSTRUMENTS

Code	Head Weight	Size
1-19790	8 oz	7.5"

GERZOG MALLET



Code	Head Weight	Size
1-19792	9 oz	7.75"

LUCAE MALLET




Code	Head Weight	Size
1-19793	7.5 oz	7.5"

MALLET WITH NYLON HEAD



Code	Description	Size
1-11312	Passer	5.5"
1-69846	Twister	6.25"


CERCLAGE WIRE



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-27160	Cuts up to 2 mm wire	6.5"

DIAMOND PIN AND WIRE CUTTER




Code	Description	Size
1-277	With Ratchet	6.0"
1-275	With Ratchet	8.5"
1-274	With Ratchet	9.5"
1-276	Without Ratchet	6.0"
1-272	Without Ratchet	9.5"

KERN BONE HOLDING FORCEPS



Code	Description	Size
1-2721	With Ratchet	13.0"
1-2723	With Ratchet	17.0"
1-2720	Without Ratchet	13.0"
1-2722	Without Ratchet	17.0"

LANE BONE HOLDING FORCEPS



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-49800	Cuts up to 3.2 mm wire, Carbide Inserts	9.5"

PIN CUTTER



Code	Description	Size
1-27158TC	Carbide Inserts	7.0"

PIN & WIRE CUTTER



Code	Size
1-27418	0.25"
1-27420	0.375"
1-27422	0.5"
1-27424	0.625"
1-27426	0.75"
1-27427	0.875"
1-27428	1.0"
1-27430	1.25"
1-27432	1.5"

HIBBS CHISEL 9.0"



ORTHOPEDIC INSTRUMENTS

Code	Size
1-27438	0.25"
1-27440	0.375"
1-27442	0.5"
1-27444	0.625"
1-27446	0.75"
1-27447	0.875"
1-27448	1.0"
1-27450	1.25"
1-27452	1.5"

HIBBS OSTEOTOME, STRAIGHT, 9.0"



Code	Description	Size
1-27488	Straight	0.25"
1-27490	Straight	0.5"
1-27492	Straight	0.75"
1-27494	Straight	1.0"
1-27499	Curved	0.25"
1-27500	Curved	0.5"
1-27502	Curved	0.75"
1-27504	Curved	1.0"

LAMBOTTE OSTEOTOME, 9.0"



Code	Width	Size
1-27820	0.25"	7.0"
1-27821	0.375"	7.5"
1-27822	0.5"	7.5"
1-27824	0.75"	8.0"
1-27826	1.0"	8.25"

KEY PERIOSTEAL ELEVATOR



ORTHOPEDIC INSTRUMENTS


Code	Width	Size
1-261410	7 mm	7.5"
1-261412	17 mm	7.5"

LANGENBECK ELEVATOR



Code	Size
1-19704	$\frac{3}{0}$
1-19706	$\frac{2}{0}$
1-19708	0
1-19710	1
1-19712	2
1-19714	3
1-19716	4
1-19718	5
1-19720	6

SPRATT BONE CURETTE 6.5"



Code	Diameter
1-157-6313	$\frac{1}{2}$ "
1-157-6316	$\frac{5}{8}$ "
1-157-6319	$\frac{3}{4}$ "
1-157-6323	$\frac{7}{8}$ "
1-157-6325	1"

GALT TREPHINES



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-27700	D/E, Curved on Flat and Sideways	10.5"
1-27702	Round and Half-Round Tapered Ends	12.0"

PUTTI BONE RASP



Code	Description	Size
1-21638	D/E, Four-Sided Flat and Convex	8.5"

FOMON RASP AND FILE



Code	Description	Size
1-27703	Half-Round Blade of 3.5 mm	8.5"

BONE FILE, BAYONET SHANK



Code	Description	Size
1-19804	4 mm x 12 mm Bite	5.75"

HARTMAN RONGEUR, STRONG CURVE



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-26190	8 mm Bite, Straight	7.0"
1-26192	8 mm Bite, Curved	7.0"

LUER RONGEUR



Code	Description	Size
1-19808	5 mm x 15 mm Bite	7.0"

BANE RONGEUR, CURVED



Code	Description	Size
1-191230	2.5 mm Bite, Straight	6.25"
1-191232	2.5 mm Bite, Curved	6.25"

LEMPERT RONGEUR



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-174800	1.3 mm Wide Jaw, Straight	5.5"
1-174801	1.3 mm Wide Jaw, Curved	5.5"

MICRO FRIEDMAN RONGEUR



Code	Size
1-271390	5.5"
1-271391	6.5"
1-271392	7.5"
1-271394	8.5"

LISTON BONE CUTTING RONGEUR



Code	Description	Size
1-25454	9 mm X 15 mm Jaw, Straight	8.75"
1-25456	9 mm X 15 mm Jaw, Curved	8.75"
1-25458	6 mm X 12 mm Jaw, Duckbill Jaw	9.0"

STILLE LUER RONGEUR, DOUBLE ACTION



ORTHOPEDIC INSTRUMENTS

Code	Description	Size
1-21626	Straight	5.5"
1-21626	Straight	7.5"
1-25386	Straight	9.0"
1-21627	Curved	5.5"
1-25385	Curved	7.5"

RUSKIN RONGEUR, DOUBLE ACTION



Code	Description	Size
1-19844	Slightly Curved	7.0"

BEYER RONGEUR, DOUBLE ACTION



Code	Description	Size
1-19724	2 mm x 10 mm Jaw	7.0"

FERRIS SMITH IVD RONGEUR, STRAIGHT



Code	Description	Size
1-271295	2 mm x 13 mm Jaw	6.0"


KLEINERT-KUTZ RONGEUR, DOUBLE ACTION



OB / GYN INSTRUMENTS


Code	Size
1-301000	10.0"

SCHROEDER UTERINE VULSELLUM FORCEPS




Code	Size
1-7620	10.5"

BOZEMAN UTERINE DRESSING FORCEPS




Code	Description	Size
1-301420	Straight	8.5"
1-301425	Curved	8.5"

SCHUBERT UTERINE BIOPSY FORCEPS



Code	Description	Size
1-28294	3.5 mm x 10 mm Bite	9.0"
1-28296	3.5 mm x 10 mm Bite	13.0"
1-28298	3.5 mm x 10 mm Bite	15.0"


YEOMAN BIOPSY FORCEPS WITH BASKET



OB / GYN INSTRUMENTS

Code	Size
1-3005	Virgin
1-3010	0.75" x 3.0"
1-3015	1.375" x 4.0"
1-3020	1.5" x 4.5"

GRAVES VAGINAL SPECULUM




Code	Size
1-3040	Virgin
1-3050	0.5" x 3.0"
1-3055	0.875" x 4.0"
1-3060	1.0" x 4.75"

PEDERSON VAGINAL SPECULUM



Code	Size
1-302350	15.5"

KIELLAND OBSTETRICAL FORCEPS



ENT INSTRUMENTS

Code	Size
1-192100	3.0"
1-19405	3.5"
2-19405	3.5"
1-19410	5.0"
2-19410	5.0"
1-19415	6.5"

ALLIGATOR FORCEPS, SERRATED JAWS



Code	Description	Size
1-19260	Size 00, Sharp	6.5"
1-19264	Size 1, Sharp	6.5"
1-19268	Size 3, Sharp	6.5"
2-19268	Size 3, Sharp	6.5"
1-19280	Size 00, Blunt	6.5"
1-19282	Size 0, Blunt	6.5"
1-19286	Size 2, Blunt	6.5"

BUCK EAR CURETTE, STRAIGHT



Code	Description	Size
1-19270	Size 00, Sharp	6.5"
1-19274	Size 1, Sharp	6.5"
1-19276	Size 3, Sharp	6.5"
1-19290	Size 00, Blunt	6.5"
1-19292	Size 0, Blunt	6.5"
2-19292	Size 0, Blunt	6.5"
1-19293	Size 1, Blunt	6.5"
1-19296	Size 2, Blunt	6.5"
2-19296	Size 2, Blunt	6.5"


BUCK EAR CURETTE, ANGLED



ENT INSTRUMENTS

Code	Description	Size
1-19314	Small Loop	6.5"
2-19314	Small Loop	6.5"
1-19316	Medium Loop	6.5"
2-19316	Medium Loop	6.5"
1-19318	Large Loop	6.5"
2-19318	Large Loop	6.5"

BILLEAU EAR LOOP



Code	Description	Size
1-181968	5 mm Blade	7.0"

FREER ELEVATOR




Code	Description	Size
1-191330	2 mm Blade, Straight	7.25"
1-191334	2 mm Blade, Curved	7.25"

LEMPERT ELEVATOR



Code	Description	Size
1-2028	2" Long Blades	6.0"
1-2030	3" Long Blades	6.0"

KILLIAN NASAL SPECULUM



ADDITIONAL INSTRUMENTS

Code	Dimensions
1-3400	8" x 5" x 2"
1-3404	12" x 7" x 2"
1-3408	16" x 9" x 2"

INSTRUMENT TRAY



Code	Dimensions
1-3400D-2	Lid for 8" x 5" x 2"
1-3404-2	Lid for 12" x 7" x 2"
1-3408-2	Lid for 16" x 9" x 2"

INSTRUMENT TRAY LIDS



Code	Size
1-1072	5.5"
2-1072	5.5"
1-1074	6.0"
2-1074	6.0"

GROOVED DIRECTOR



ADDITIONAL INSTRUMENTS

Code	Qty.
450205	50/ Carton

WEC PREP BLADES



Code	Description	Size
381-201	Babcock Suture Wire	18 Gauge
381-202	Babcock Suture Wire	20 Gauge
381-228	Babcock Suture Wire	22 Gauge
381-269	Babcock Suture Wire	26 Gauge
381-285	Babcock Suture Wire	28 Gauge
381-301	Babcock Suture Wire	30 Gauge

STAINLESS STEEL ORTHOPEDIC WIRE



Code	Color	Size
T-01-01	White	0.125"
T-01-02	Blue	0.125"
T-01-03	Green	0.125"
T-01-04	Red	0.125"
T-01-05	Yellow	0.125"
T-01-07	Orange	0.125"
T-01-09	Black	0.125"

INSTRUMENT ID TAPE 250"



PILLING® & KMEDIC® INSTRUMENTS

Code
465160

Size
5.625"

THUMB FORCEPS 2 mm TIP



Code
057150

Size
7.0" (4.0" Stem)

ALLIGATOR FORCEPS



Code
KM47870

Description
4 mm x 10 mm Bite

Size
7.0"

FERRIS SMITH RONGEUR, STRAIGHT



PILLING® & KMEDIC® INSTRUMENTS

Code	Size
181816	5.0"

PLASTI-STAT MOSQUITO FORCEPS CURVED



Code	Description	Size
510440	Carbide Inserts	5.5"
510445	Carbide Inserts	6.5"

OLSEN-HEGAR NEEDLE HOLDER



Code	Size
KM71111	9.0" x 0.035"
KM71112	9.0" x 0.045"
KM71113	9.0" x 0.062"

K-WIRE DOUBLE END TROCAR (6/PK)



KMEDIC® SUPER-CUT SCISSORS

- **The industry's sharpest scissors**, saving time and eliminating frustration.
- Cut through soft and fibrous tissue without effort.
- Less tearing during dissection/cutting means faster healing and less scarring.
- **Retains sharpness much longer than ordinary scissors.**
- Less hand fatigue for the surgeon and less trauma for the patient.
- KMedic Certified Instruments are superior instruments providing lasting value.

Code	Size
KM35202	5.5"

SUPER-CUT HI-LEVEL BANDAGE SCISSORS



Code	Description	Size
KM35210		5.5"
KM35212		7.25"
KM35214	With Large Ring	8.0"

SUPER-CUT LISTER BANDAGE SCISSORS



Code	Description	Size
KM35220	Straight	5.5"
KM35250	Curved	5.5"

SUPER-CUT OPERATING SCISSORS SHARP/BLUNT



KMEDIC SUPER-CUT SCISSORS

Code	Description	Size
KM35280	Straight	5.5"
KM35292	Straight	6.75"
KM35316	Straight	9.0"
KM35286	Curved	5.5"
KM35298	Curved	6.75"
KM35322	Curved	9.0"
KM35340	Curved	11.0"

SUPER-CUT MAYO SCISSORS



Code	Description	Size
KM35364	Straight, Baby Metz	4.5"
KM35380	Straight	6.0"
KM35392	Straight	7.0"
KM35370	Curved, Baby Metz	4.5"
KM35382	Curved	6.0"
KM35394	Curved	7.0"
KM35396	Curved, Delicate	7.0"
KM35406	Curved	8.0"
KM35418	Curved	9.0"
KM35430	Curved	11.0"

SUPER-CUT METZENBAUM SCISSORS



Code	Description	Size
KM35450	Straight, Delicate	3.5"
KM35454	Straight	4.5"
KM35452	Curved	3.5"
KM35460	Curved, Delicate	4.5"

SUPER-CUT IRIS SCISSORS



KMEDIC SUPER-CUT SCISSORS

Code	Description	Size
KM35484	Straight	4.5"
KM35490	Curved	4.5"
KM35494	Curved	5.25"

SUPER-CUT STEVENS SCISSORS, BLUNT POINTS



Code	Description	Size
KM35496	Curved	6.25"

SUPER-CUT JAMISON SCISSORS



SURGICAL INSTRUMENT PACKS

Select surgical instrument packs to simplify your purchase and save you money.
Available in two grades.

SPAY PACK

Includes the following 18 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder	5.5"	1	1-815	2-815
Brown Adson Tissue Forceps	4.75"	1	1-6124	2-6124
Tissue Forceps, 1 x 2 Teeth	5.5"	1	1-644	2-644
Backhaus Towel Clamp	3.5"	4	1-7504	2-7504
Mayo Scissors, Straight	6.75"	1	1-5124	2-5124
Metzenbaum Scissors, Curved	7.0"	1	1-5182	2-5182
Halsted Mosquito Forceps, Curved	5.0"	2	1-74	2-74
Kelly Forceps, Curved	5.5"	2	1-738	2-738
Rochester Carmalt Forceps, Curved	6.25"	2	1-7172	2-7172
Allis Tissue Forceps, 4 x 5 Teeth	6.0"	1	1-166	2-166
Knife Handle #3	5.0"	1	1-47	2-47
Spay Hook	8.0"	1	1-28131	2-28131
To order the complete pack, use the following list number:			SP-1	SP-2

GENERAL SURGERY PACK

Includes the following 19 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder, T/C	5.5"	1	1-815TC	2-815TC
Brown Adson Tissue Forceps	4.75"	1	1-6124	2-6124
Backhaus Towel Clamp	5.5"	4	1-7506	2-7506
Mayo Scissors, Straight	6.75"	1	1-5124	2-5124
Metzenbaum Scissors, Curved	7.0"	1	1-5182	2-5182
Halsted Mosquito Forceps, Curved	5.0"	3	1-74	2-74
Kelly Forceps, Straight	5.5"	2	1-736	2-736
Rochester Carmalt Forceps, Curved	6.25"	4	1-7172	2-7172
Knife Handle #3	5.0"	1	1-47	2-47
Spay Hook	8.0"	1	1-28131	2-28131
To order the complete pack, use the following list number:			GSP-1	GSP-2

SURGICAL INSTRUMENT PACKS

CANINE SPAY PACK

Includes the following 19 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder	6.5"	1	1-816	2-816
Brown Adson Tissue Forceps	4.75"	1	1-6124	2-6124
Backhaus Towel Clamp	3.5"	4	1-7504	2-7504
Mayo Scissors, Straight	6.75"	1	1-5124	2-5124
Halsted Mosquito Forceps, Straight	5.0"	2	1-72	2-72
Halsted Mosquito Forceps, Curved	5.0"	2	1-74	2-74
Kelly Forceps, Straight	5.5"	2	1-736	2-736
Kelly Forceps, Curved	5.5"	2	1-738	2-738
Allis Tissue Forceps, 4 x 5 Teeth	6.0"	1	1-166	2-166
Knife Handle #3	5.0"	2	1-47	2-47
Spay Hook	8.0"	1	1-28131	2-28131
To order the complete pack, use the following list number:			CSP-1	CSP-2

FELINE SPAY PACK

Includes the following 17 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder	5.5"	1	1-815	2-815
Brown Adson Tissue Forceps	4.75"	1	1-6124	2-6124
Backhaus Towel Clamp	3.5"	4	1-7504	2-7504
Mayo Scissors, Straight	6.75"	1	1-5124	2-5124
Halsted Mosquito Forceps, Straight	5.0"	2	1-72	2-72
Halsted Mosquito Forceps, Curved	5.0"	2	1-74	2-74
Kelly Forceps, Straight	5.5"	2	1-736	2-736
Kelly Forceps, Curved	5.5"	2	1-738	2-738
Knife Handle #3	5.0"	1	1-47	2-47
Spay Hook	8.0"	1	1-28131	2-28131
To order the complete pack, use the following list number:			FSP-1	FSP-2

SURGICAL INSTRUMENT PACKS

STARTER OPHTHALMIC PACK

Includes the following 4 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Derf Needle Holder	4.75"	1	1-84	2-84
Bishop Harmon Forceps, 1 x 2 Teeth	3.4"	1	1-18867	2-18867
Barraquer Wire Speculum (Large Blades)	1.75"	1	1-1835	2-1835
Stevens Tenotomy Scissors, Curved	4.2"	1	1-181466	2-181466
To order the complete pack, use the following list number:			SOP-1	SOP-2

ADVANCED OPHTHALMIC PACK

Includes the following 7 instruments and 1 tray:

Description	Size	Quantity	Grade 1	Grade 2
Bishop Harmon Forceps, 1 x 2 Teeth	3.4"	1	1-18867	2-18867
Desmarres Chalazion Forceps	3.5"	1	1-181202	2-181202
Barraquer Wire Speculum (Large Blades)	1.75"	1	1-1835	2-1835
Stevens Tenotomy Scissors, Curved	4.2"	1	1-181466	2-181466
Castroviejo Needle Holder With Lock, Straight	5.5"	1	1-181828	2-181828
Handle For Miniature Edged Blade	—	1	1-1310	**
#64 Miniature Edged Blade	—	12/Box	1-176400	**
SteriBest® Instrument Tray With Mat	6.5" x 4.0"	1	CP634	**
To order the complete pack, use the following list number:			AOP-1	AOP-2

KNEE PACK

Includes the following 8 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Mini-Hohmann Retractor, 6mm Blade	6.25"	1	1-27213	**
Mini-Hohmann Retractor, 8mm Blade	6.25"	1	1-27214	**
Hohmann Retractor, 18mm Blade	9.25"	1	1-27216	**
Halsted Mosquito Forceps, 1x2 Teeth, Straight	5.0"	1	1-714	2-714
Senn Retractor, Sharp	6.25"	2	1-1174	**
Gelpi Retractor	3.5"	1	1-11364	2-11364
Gelpi Retractor	5.5"	1	1-11365	2-11365
To order the complete pack, use the following list number:			KNEE-1	KNEE-2

SURGICAL INSTRUMENT PACKS

EAR PACK

Includes the following 11 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Spratt Bone Curette, Size 0	6.5"	1	1-19708	**
Olsen-Hegar Needle Holder	5.5"	1	1-815	2-815
Brown Adson Tissue Forceps	4.75"	1	1-6124	2-6124
Metzenbaum Scissors, Curved	7.0"	1	1-5182	2-5182
Halsted Mosquito Forceps, Curved	5.0"	2	1-74	2-74
Dean Rongeur, Curved	6.0"	1	1-22480	**
Weitlaner Retractor, Sharp, 2 x 3 Teeth	4.5"	1	1-11600	**
Weitlaner Retractor, Sharp, 3 x 4 Teeth	5.5"	1	1-11608	**
Knife Handle #3	5.0"	1	1-47	2-47
Senn Retractor, Sharp	6.25"	1	1-1174	**
To order the complete pack, use the following list number:			EAR-1	EAR-2

EQUINE SURGERY PACK

Includes the following 19 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder	6.5"	2	1-816	2-816
Dressing Forceps, Serrated Tips	6.0"	1	1-610	2-610
Brown Adson Tissue Forceps	4.75"	2	1-6124	2-6124
Adson Dressing Forceps, Serrated	4.75"	1	1-6118	2-6118
Lorna Towel Clamp	5.25"	4	1-7515	2-7515
Mayo Scissors, Straight	6.75"	2	1-5124	2-5124
Metzenbaum Scissors, Straight	7.0"	1	1-5181	2-5181
Halsted Mosquito Forceps, Straight	5.0"	2	1-72	2-72
Rochester Carmalt Forceps, Straight	8.0"	2	1-7180	2-7180
Ragnell Retractor	6.0"	1	1-1173	2-1173
Knife Handle #3	5.0"	1	1-47	2-47
To order the complete pack, use the following list number:			ESP-1	ESP-2

SURGICAL INSTRUMENT PACKS

EQUINE LONG INSTRUMENT PACK

Includes the following 8 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olson-Hegar Needle Holder	7.5"	1	1-817	2-817
Stratte Needle Holder, T/C	9.0"	1	1-8102TC	**
Dressing Forceps, Serrated Tips	10.0"	1	1-616	2-616
Metzenbaum Scissors, Straight	11.0"	1	1-5190	**
Allis Tissue Forceps, 5 x 6 Teeth	10.0"	1	1-1612	**
Rochester Pean Forceps, Curved	9.0"	1	1-7144	**
Knife Handle #3L	8.375"	1	1-410	2-410
Knife Handle #4L	8.375"	1	1-411	2-411
To order the complete pack, use the following list number:			ELP-1	ELP-2

EQUINE THROAT PACK

Includes the following 23 instruments:

Description	Size	Quantity	Grade 1	Grade 2
Olsen-Hegar Needle Holder	7.5"	1	1-817	2-817
Dressing Forceps, Serrated Tips	10.0"	2	1-616	2-616
Metzenbaum Scissors, Straight	11.0"	2	1-5190	**
Mayo Scissors, Curved	10.0"	1	1-5134	**
Sims Scissors, Curved, S/B	8.0"	1	1-5222	**
Crile Forceps, Curved	6.5"	2	1-748	2-748
Allis Tissue Forceps, 5 x 6 Teeth	10.0"	3	1-1612	**
Rochester Pean Forceps, Straight	9.0"	3	1-7126	**
Rochester Pean Forceps, Curved	9.0"	2	1-7144	**
Babcock Intestinal Forceps, 15mm Jaws	9.25"	1	1-1648	**
Foerster Sponge Forceps, Straight	9.5"	1	1-7600	2-7600
Volkman Curette, D/E, Oval Cups	5.5"	1	1-27920	**
Gelpi Retractor	6.75"	1	1-11366	2-11366
Rigby Retractor With Grip Lock	6.75"	1	1-11374	**
Knife Handle #3L	8.375"	1	1-410	2-410
To order the complete pack, use the following list number:			ETP-1	ETP-2

SURGICAL INSTRUMENT PACKS

MICROSURGICAL INSTRUMENT PACK

Includes the following 9 instruments and 1 tray:

Description	Size	Quantity	Grade 1	Grade 2
Microsurgery Needle Holder, Pencil Grip, w/o Lock, Straight, T/C	6.5"	2	1-171026TC	**
Tewes-Barraquer Micro Dissecting Scissors, Pencil Grip, Straight	6.75"	1	1-394441	**
Micro Tying Forceps, 0.5mm Wide Tip, Round Handle, Straight	7.125"	1	357660	**
Micro Tying Forceps, 0.5mm Wide Tip, Round Handle, Curved	7.125"	1	357661	**
Micro Ring Tip Forceps, 0.5mm Wide Tip, Pencil Grip, Straight	7.0"	1	1-357647	**
Knife Handle #3	5.0"	1	1-47	2-47
Strabismus Scissors, Straight	4.0"	1	1-5312	**
Adson Dressing Forceps, Straight	4.75"	1	1-6118	2-6118
Alm Retractor	2.5"	1	1-1112	2-1112
SteriBest® Instrument Tray With Mat	10.0 x 6.0"	1	CP10385	**
To order the complete pack, use the following list number:			MIP-1	MIP-2

If you have a specific instrument pack you would like to have developed, please contact RICA's office at (800) 889-3218 for assistance.

INSTRUMENT TRAYS

With Pin Mats or Silicone Bar Holders

SteriBest® Instrument Trays are designed to protect delicate instruments during sterilization, storage, and transportation. Available in various sizes and styles, they feature silicone mats or holders for added protection. Made in the USA from durable, heat- and chemical-resistant Ultem®, these trays are compatible with steam and gas-plasma sterilization. Raised nodes and enhanced ventilation promote faster drying, while transparent, slotted lids allow easy content identification and improve steam penetration to prevent “wet pack” issues.


- Advanced design allows for faster sterilization and drying.
- Designed to protect delicate instruments during sterilization and storage.
- Made from the polymer Ultem® to ensure strength and long lasting durability.
- Transparent lids for easy identification of contents.
- Trays are stackable for storage.





ADVANCED DESIGN FEATURES:

Slotted lids and bases

facilitate circulation and more rapid drying.



Nodes on bases permit better steam circulation under the mats.



INDEPENDENT TEST LAB VALIDATIONS:

DRY TIME STUDY*

Independent lab dry time results validated 25 minutes on all trays, giving you guidance to set your parameters to minimize wet pack.

*Based on standards used in the studies.

STERILIZATION VALIDATION STUDY*

Independent testing showed 100% kill both above and below the mats at 2 minutes, giving you confidence that your instrumentation will be sterile.

Ultem® is a registered trademark of SABIC Innovative Plastics, Inc. #85198347

INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 6" LENGTH

Code	Dimensions
CP614	6.0" x 2.5" x .75" 152mm x 64mm x 19mm



TRAY WITH BASE, LID & BARS



CP614-DK-1



CP614-DK-2




CP614-DK-3

Code	Description	Dimensions
CP614-DK-1	2 Bars with 1-Slot Knife Holder	6.0 in. x 2.5 in. x .75 in. 5.5mm bar slot
CP614-DK-2	2 Bars with 2-Slot Knife Holder	6.0 in. x 2.5 in. x .75 in. 5.5mm bar slots
CP614-DK-3	2 Bars with 3-Slot Knife Holder	6.0 in. x 2.5 in. x .75 in. 5.5mm bar slots

INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 7.5" LENGTH

Code	Dimensions
CP778	7.5 in. x 2.5 in. x .75 in. 190.5mm x 63.5mm x 19mm



TRAY WITH BASE, LID & BARS



CP778-DK-1



CP778-DK-2




CP778-DK-3

Code	Description	Dimensions
CP778-DK-1	2 Bars with 1-Slot Knife Holder	7.5 in. x 2.5 in. x .75 in. 5.5mm bar slot
CP778-DK-2	2 Bars with 2-Slot Knife Holder	7.5 in. x 2.5 in. x .75 in. 5.5mm bar slots
CP778-DK-3	2 Bars with 3-Slot Knife Holder	7.5 in. x 2.5 in. x .75 in. 5.5mm bar slots

INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 6.5" Length

Code	Dimensions
CP634	6.5 in. x 4.0 in. x .75 in. 165mm x 102mm x 19mm




TRAY WITH BASE, LID & BARS

Code	Description	Dimensions
CP634-DK-5	2 Bars with 5-Slot Knife Holder	6.5 in. x 4.0 in. x .75 in. 5.5mm bar slot



TRAY WITH BASE, LID, MAT, INSERT TRAY & INSERT MAT 10" LENGTH


Code	Dimensions
CP1038D2	10.0 in. x 6.0 in. x 1.5 in. 254mm x 152mm x 38mm



INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 10" LENGTH

Code	Dimensions
CP1038D1	10.0 in. x 6.0 in. x 1.5 in. 254mm x 152mm x 38mm



TRAY WITH BASE, LID & BARS



CP1038D1-P-1S/L




CP1038D1-P-3S/L

Code	Description	Dimensions
CP1038D1-P-1S/L	2 Bars with 1-Slot Phaco Holder	10.0 in. x 6.0 in. x 1.5 in. Bar slots: 16mm; 8mm
CP1038D1-P-3S/L	2 Bars with 3-Slot Phaco Holder	10.0 in. x 6.0 in. x 1.5 in. Bar slots: 7.5mm; 15.5mm, 7.5mm and 12.5mm

TRAY WITH BASE, LID & MAT 10" LENGTH


Code	Dimensions
CP1038S	10.0 in. x 6.0 in. x .75 in. 254mm x 152mm x 19mm



INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 7.5" LENGTH

Code	Dimensions
CP471D1	7.5 in. x 4.0 in. x 1.5 in. 190.5mm x 102mm x 38mm



TRAY WITH BASE, LID & BARS



CP471D1-P-1S/L




CP471D1-P-3S/L

Code	Description	Dimensions
CP471D1-P-1S/L	2 Bars with 1-Slot Phaco Holder	7.5 in. x 4.0 in. x 1.5 in. Bar slots: 16mm; 8mm
CP471D1-P-3S/L	2 Bars with 3-Slot Phaco Holder	7.5 in. x 4.0 in. x 1.5 in. Bar slots: 7.5mm; 15.5mm, 7.5mm and 12.5mm

INSTRUMENT TRAYS

TRAY WITH BASE, LID & MAT 7.5" LENGTH

Code	Dimensions
CP471S	7.5 in. x 4.0 in. x .75 in. 190.5mm x 102mm x 19mm



TRAY WITH BASE, LID & BARS



CP471S-DK-2



CP471S-DK-3



CP471S-DK-4



CP471S-DK-5

Code	Description	Dimensions
CP471S-DK-2	2 Bars with 2-Slot Knife Holder	7.5 in. x 4.0 in. x .75 in. 5.5mm bar slots
CP471S-DK-3	2 Bars with 3-Slot Knife Holder	7.5 in. x 4.0 in. x .75 in. 5.5mm bar slots
CP471S-DK-4	2 Bars with 4-Slot Knife Holder	7.5 in. x 4.0 in. x .75 in. 5.5mm bar slots
CP471S-DK-5	2 Bars with 5-Slot Knife Holder	7.5 in. x 4.0 in. x .75 in. 5.5mm bar slots

INSTRUMENT ENZYMATIC CLEANING SOLUTION



- Penetrates hard to reach areas to help dissolve, loosen and flush blood, protein soils, and tissue.
- Eliminates hand scrubbing.
- Non-corrosive with a neutral PH.
- Fresh clean scent.
- Can be used as an instrument pre-soak.
- Excellent concentrate to use in ultrasonic cleaners.
- Fast acting, specialized enzymatic cleaner.
- Manufactured in a cGMP compliant facility for batch-to-batch consistency.
- Economical and easy to use.

Code	Description
REC128	1 Gallon Bottle With Pump



AUTOCLAVE GLOVES

- 100% cotton, terry cloth.
- Absorbs moisture.
- Heat resistant, up to 250°F.
- Launderable.
- Natural color with black contrast edge.
- Available in 15.5" and 17.0" lengths.
- One size fits all hands.



Code

ACG15

ACG17

Description

15.5" length to cover the wrist

17.0" length for protection up to the elbow

CLEANING BRUSHES



Micro Channel



Small Diameter



Medium/Large Diameter



Diamond Shape



80.TB-2100



80.TN-2121

Code	Description	Bead	Brush Diameter	Length	Pkg Qty
80.M-08-005	Micro channel, pull through	-	1.5 FR / 0.02" / .50 mm	8"	3
80.M-12-020	Micro channel, pull through	-	6 FR / 0.08" / 2.0 mm	12"	3
80.FR-12-03	Small diameter, pull through	Purple	3 FR / 0.04" / 1.0 mm	12"	3
80.FR-12-04	Small diameter, pull through	Aqua	4 FR / 0.06" / 1.4 mm	12"	3
80.FR-12-05	Small diameter, pull through	Red	5 FR / 0.07" / 1.7 mm	12"	3
80.FR-12-07	Small diameter, pull through	Blue	7 FR / 0.09" / 2.3 mm	12"	3
80.B-16-030	Small diameter, pull through	Yellow	0.12" / 3.0 mm	16"	3
80.B-16-040	Small diameter, pull through	Orange	0.16" / 4.0 mm	16"	3
80.B-16-050	Medium/large diameter, pull through	Green	0.20" / 5.0 mm	16"	3
80.B-20-050	Medium/large diameter, pull through	Green	0.20" / 5.0 mm	20"	3
80.B-16-120	Medium/large diameter, pull through	Brown	0.47" / 12.0 mm	16"	3
80.B-16-100	Medium/large diameter, pull through	White	0.39" / 10.0 mm	16"	3
80.D-10-159	Diamond shape, pull through	Green	0.22 x 0.63 x 0.22" (5.5 x 15.9 x 5.5 mm)	10"	3
80.D-16-103	Diamond shape, pull through	White	0.25 x 0.41 x 0.25" (6.4 x 10.3 x 6.4 mm)	16"	3
80.TB-2100	Red toothbrush, brass bristles	-	1/4"	7"	10
80.TN-2121	Dark taupe toothbrush, nylon bristles, double-ended	-	1/4"	7"	10

INSTRUMENT CLEANING

STAIN REMOVER

Code	Description	Pkg Qty
80.KSLSR-8501	Instrument Stain Remover - 3 oz jar	1
80.KSLSR-8512	Instrument Stain Remover - 3 oz jar	12



ADHESIVE REMOVER

Code	Description	Pkg Qty
80.T-SAMP-RMV	Tape-B-Gone Adhesive Tape Remover, 0.5 oz bottle	1
80.T-08RMV	Tape-B-Gone Adhesive Tape Remover, 8 oz bottle	1
80.T-32RMV	Tape-B-Gone Adhesive Tape Remover, 8 oz bottle	4

CLEANING PADS

Code	Description	Colour	Dimensions	Pkg Qty
80.3FRIN556120	Rough	Brown	6 x 4"	1
80.3FRIN556125	Fine	Gray	6 x 4"	12



FACIAL PROTECTION



FULL FACE SHIELD

Solid or vented foam for breathability.
Soft elastic or adjustable Velcro straps, secured with an eyelet for added safety.



PROTECTIVE EYEWEAR

Non-perforated with self-closing holes.
Top-down fluid protection ridge on frame.
Wrap-around design prevents slippage.



80.11000-507



80.11000-112



80.11000-410

Code	Description	Pkg Qty
80.11000-000	Splash Mask w/ Velcro Strap, Vented Foam	25
80.11000-002	Splash Mask w/ Velcro Strap, Solid Foam	25
80.11000-003	Splash Mask w/ Elastic Strap, Solid Foam	25
80.11000-500	Cool ShieldZ Pre-Assembled, 10/bag-10 bags/cs	10
80.11000-506	Cool ShieldZ, Dispenser Box 20/box - 6 boxes/cs	25
80.11000-507	Wall Mount Bracket	25
80.11000-112	Solid Side Shield, Pair	25
80.11000-400	Surgical, Tie Mask w/ Attached Shield	25
80.3FSZB790111	Anti fog cloth	20

MARKERS



Code	Description	Pkg Qty
80.K1492-20	Viscot Dual Tip - Gentian Ink Regular Tip & Permanent Marker Ink with Fine Tip	20
80.K4387-20	Viscot Surgical Marker, Regular Tip	20
80.K1811-20	Viscot Waterproof Permanent Marking Pen Fine Tip (bulk)	20

WECK® VISISTAT® SKIN STAPLING SYSTEM Visistat®

- Easy-to-use device with quality performance for external skin closure.
- Loaded with 35 regular or wide medical grade stainless steel staples.
- Features a see-through window for better visibility.
- Patented non-stop feed mechanism for easy release and reduced jamming.

The Visistat Skin Stapler is perfectly supplemented by our Weckstat™ Skin Staple Extractor, providing both comfort and simplicity during staple removal.



528135



528235



525980

Code	Description	Staple Size	Quantity
528135	Visistat Skin Stapler - 35R	Regular - 5.7mm Wide by 3.9mm Tall	6 Staplers/Box
528235	Visistat Skin Stapler - 35W	Wide - 6.5mm Wide by 4.7mm Tall	6 Staplers/Box
525980	Weckstat™ Skin Staple Extractor		12 Staple Extractors/Box

HORIZON™ LIGATION CLIP SYSTEM

HORIZON®

Weck's Horizon System provides veterinary surgeons all of the benefits and time-savings associated with the Hemoclip System, with these additional benefits:

- Easier to load clip with a near frictionless one-hand load.
- Adhesive backing on the cartridge allows for placement in a convenient location.
- Additional clip size...Microclip...for accurate and secure ligation of very small blood vessels.



HORIZON APPLIER

Code	Description	Color	Length
537061	Microclip, Curved (20°)	White	6.0"
537081	Microclip, Curved	White	8.0"
137061	Small, Curved	Yellow	6.0"
137081	Small, Curved	Yellow	8.0"
237061	Medium, Curved	Blue	6.0"
237081	Medium, Curved	Blue	8.0"
237117	Medium, Right Angled (70°)	Blue	11.0"
337081	Medium/Large, Curved	Green	8.0"
337085	Medium/Large, Angled (55°)	Green	8.0"
337117	Medium/Large, Right Angled	Green	11.0"
437081	Large, Curved	Orange	8.0"
437085	Large, Angled	Orange	8.0"
437117	Large, Right Angled	Orange	11.0"

HORIZON™ LIGATION CLIP SYSTEM

HORIZON®



HORIZON CARTRIDGE

Code	Description	Color	Clips/ Cartridge	Cartridges/Box
005200	Microclip	White	6	30
001200	Small	Yellow	6	30
001204	Small	Yellow	24	25
002200	Medium	Blue	6	30
002204	Medium	Blue	24	25
003200	Medium/Large	Green	6	20
003204	Medium/Large	Green	24	15
004200	Large	Orange	6	20
004204	Large	Orange	24	15

BFW™ VETERINARY HEADLIGHTS



KITSILANO™ CORDLESS HEADLIGHT SYSTEM

- Up to ~40,000 LUX.
- Eight adjustable brightness levels plus a Turbo mode.
- Wireless, rechargeable.
- Provides 5 hours of continuous light at full intensity.
- Ultra-lightweight design.

CHROMAPRO™ LED LIGHT SOURCE

- Up to ~70,000 LUX.
- 50-watt LED illumination.
- LED expected lifespan: 30,000 hours.



BFW9030M

DOVER™ & MAUI DOVER™ HEADLIGHT SYSTEMS

- Up to ~105,000 LUX.
- Wireless, rechargeable, surgical quality.
- Two battery packs with a recharging station.

Code

BFWP300

BFW5500

BFW9030

BFW9030B

BFW9030M

Description

Kitsilano™ Cordless Headlight System

ChromaPRO™ 50 Watt LED Light Source

Dover™ Headlight System, Classic Headband, Single Battery System, Holster

Dover™ Headlight System, Classic Headband, Dual Battery System, Holster

Dover™ Headlight System, Maui Headband, Dual Battery System, Cordless

FIREFLY GLOBAL™ WIRELESS HD VIDEO OTOSCOPE

- Wireless HD Imaging – View and capture images on your device.
- Clear Otoscopy – High-quality ear examinations.
- Versatile Use – Great for dental, skin, nasal, and rectal exams.
- 35x Magnification – Sharp, detailed images.



DE571



ES140



ES108

Code

DE571

ES140

ES108

Description

Wireless HD Veterinary Video Otoscope


Flexible Forceps

XL Specula Set - Includes 6 metal speculas of various lengths

MEDIFLEX[®] DYNATRAC[®] RETRACTOR SYSTEM

- Ensures optimal surgical exposure and access.
- Supports superficial and deep retraction with Surgical Stays, rakes, and blades.
- Eliminates the need for assistant-held retraction.
- Uses fewer Stays per procedure, reducing costs.
- Repposable System balances cost savings and sustainability.

Mediflex[®]
SURGICAL PRODUCTS



MEDIFLEX® DYNATRAC® RETRACTOR SYSTEM

REUSABLE

DynaTrac® Frames

Mediflex®
SURGICAL PRODUCTS

HORSESHOE



SQUARE



FIGURE 8



DISPOSABLE

DynaTrac® Silicone Surgical Stays



**FULL RANGE OF REUSABLE STAINLESS STEEL
RETRACTOR BLADES WITH TILTING RATCHES
ARE ALSO AVAILABLE.**

Code	Description	Pkg Qty
78000-HP	Dynatrac Retractor Frame, hex post, horseshoe, 8.4" x 7.8"	1
78001-HP	Dynatrac Retractor Frame, hex post, square, 7.4" x 7.4"	1
78002-HP	Dynatrac Retractor Frame, 2 hex post, figure 8, 15.3" x 8.9"	1
73000	Strongarm, single, front mount, QC fitting and rail clamp (optional)	1
4005B	Surgical Stays, 5mm blunt hook	4
4005S	Surgical Stays, 5mm sharp hook	4
4012B	Surgical Stays, 12mm blunt hook	4

TECNIMED VISIOFOCUS VET THERMOMETER

- No-touch infrared temperature readings.
- Multi-site measurement from the eye, auricle, gums, or rectal area.
- Hygienic & precise for veterinary use.
- Syncs with VisioVet app for any animal.
- Adapts to environmental conditions.
- Suitable for all animals.



Customized Correction Settings To Suit Different Animals



Code


06610

Description

Infrared non-contact thermometer

PAW PROTECTOR® CAGE RACKS

- Keep animals dry and clean by keeping animals 1 ½" from the cold cage floor.
- Strong, safe, durable and easy to clean.
- Prevents injured paws as the pet's toe cannot become caught, twisted or injured in the cage rack.
- Constructed of a superior grade of strong, durable polypropylene plastic.
- Sits flat on the cage floor for a stable surface.
- **10 year warranty** against chewing and breakage.
- Available in several sizes and will fit snugly in Shor-Line®, Suburban® & Midmark® brand cages.



Code	Cage Dimensions	Rack Dimensions	Weight
1018/Rev2	Outside 18" x 28", Inside 15.5" x 27.5"	15.25" x 27.25"	2.7 lbs.
1024/Rev2	Outside 24" x 28", Inside: 21.5" x 27.5"	21.25" x 27.25"	3.8 lbs.
1030/Rev 2	Outside 30" x 28", Inside 27.5" x 27.5"	27.25" x 27.25"	4.8 lbs.
1036/Rev 2	Outside 36" x 28", Inside 33.5" x 27.5"	33.25" x 27.25"	5.9 lbs.
1048/Rev 2*	Outside 48" x 28", Inside 45.5" x 27.5"	45.25" x 27.25"	7.9 lbs.
1060/Rev 2*	Outside 60" x 28", Inside 57.5" x 27.5"	57.25" x 27.25"	10.2 lbs
1072/Rev 2**	Outside 72" x 28", Inside 69.5" x 27.5"	69.25" x 27.25"	12.4 lbs

*All of the racks identified with this symbol are actually two racks that sit side-by-side each other in the individual cage. Each rack of the pair is identified with laser etching to make it easy to know which pieces go together. For example, the 1048/Rev2 size rack is identified as part 1048A (25.61" x 27.25" size) and 1048B (19.64" x 27.25" size).

**The 1072/Rev2 rack is actually three racks that sit side-by-side each other in the individual cage. Each section of the rack is identified as 1072A, 1072B and 1072C, so you know the parts to this size rack.

HOW TO ORDER

To place an order to purchase any RICA Animal Health instruments and products or to inquire about the instruments available and current pricing, please contact your local sales representative, visit RICA's website at www.RICAAnimalhealth.com, or contact RICA Animal Health's office located in Schiller Park:

9207 Ivanhoe Street
Schiller Park, Illinois 60176
1-800-889-3218
Fax: 847-671-5690
Email: sales@RICAAnimalhealth.com

PAYMENTS

RICA Animal Health accepts Mastercard, VISA, American Express, and Discover credit cards for payment of invoices. Payment terms are 30 days net from the date of invoice. An interest charge of 1-1/2% per month may be assessed on past due invoices.

Make payment payable to RICA Animal Health and send payment to:

RICA Animal Health
Accounts Receivable
9207 Ivanhoe Street
Schiller Park, IL 60176

RETURN GOODS POLICY

All returns must be pre-approved by RICA Animal Health, which can be done by contacting us at 1-800-889-3218. Authorized returns must be made within thirty (30) days after receipt of the order. Instruments must be unused and in the original packaging. There will be a 20% restocking fee on returns over 30 days but before 90 days, a 30% restocking fee on returns over 90 days but before 120 days. No returns will be accepted 120 days after receipt of shipment. Sterile products, which have been opened, returned out of the original package or returned after the sterile expiration date are not returnable. Return goods must be shipped freight prepaid to RICA Animal Health, 9207 Ivanhoe Street, Schiller Park, Illinois 60176 and be accompanied by a copy of the original packing list and/or invoice. All items are subject to RICA's inspection and acceptance.

The RICA office is staffed from 8:00 a.m. to 4:30 p.m. CST. Orders are shipped via United Parcel Service ground service. The actual shipping charges will be added to the invoiced amount. If you require expedited delivery, please let us know when you place your order.

STERIBEST[®] PRODUCT LIMITED WARRANTY

All SteriBest[®] products are covered by RICA Animal Health's limited warranty. RICA Animal Health (RICA) warrants this product against defects in material or workmanship. RICA will repair or replace the defective product at its option if returned within the warranty period to our service facility. This warranty is in lieu of all other warranties, expressed or implied, including, but not limited to any implied warranty of merchantability or fitness for a particular purpose. The remedies provided herein are the buyer's sole and exclusive remedies. RICA shall not be liable for any direct, indirect, special, incidental, or consequential damages, whether based on contract, tort, or any legal theory. In the United States, some states do not allow exclusion or limitation of incidental or consequential damages, so the limitations above may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. This section does not limit liability for bodily injury of a person.

RICA's entire liability for any defective product shall in no event exceed the purchase price for the defective product. This limitation applies even if RICA cannot or does not repair or replace any defective product and your exclusive remedy fails of its essential purpose. See the RICA website for full details of limited warranty. Contact RICA Animal Health for more information.

INDEX

1-6826S	6	1-5179TC	9	1-514	12	1-9122	15
1-5120	7	1-5181	9	1-516	12	1-5510	16
2-5120	7	2-5181	9	2-516	12	1-5512	16
1-5120TC	7	1-5181TC	9	1-516TC	12	2-5512	16
1-5136R	7	1-5186	9	1-518	12	1-5514	16
1-5124	7	2-5186	9	1-52	12	2-5514	16
2-5124	7	1-5186TC	9	1-54	12	1-5516	16
1-5124TC	7	1-5190	9	2-54	12	2-5516	16
1-103-5017L	7	1-5190TC	9	1-56	12	1-5516LH	16
1-132-860L	7	1-5178	10	2-56	12	1-5550	16
1-5128	7	1-5180	10	1-56TC	12	2-5550	16
1-5128TC	7	2-5180	10	1-522	12	1-5998	16
1-5132	7	1-5180TC	10	1-524	12	1-51000	16
1-5132TC	7	1-5180ATC	10	2-524	12	1-5560	17
1-5136	7	1-5286	10	1-526	12	2-5560	17
1-5136TC	7	1-104-3514L	10	2-526	12	1-5560	17
1-5122	7	1-132-711L	10	1-526TC	12	1-5578	17
2-5122	7	1-5182	10	1-542	13	1-5564	17
1-5122TC	7	2-5182	10	1-544	13	2-5564	17
1-5118	7	1-5182TC	10	2-544	13	1-64	18
1-5126	7	1-132-750L	10	1-546	13	1-66	18
2-5126	7	1-132-751L	10	2-546	13	2-66	18
1-5126TC	7	1-5184	10	1-546TC	13	1-610	18
1-103-5117L	7	1-5184TC	10	1-549	13	2-610	18
1-132-820L	7	1-5188	10	1-532	13	1-614	18
1-5130	7	2-5188	10	2-532	13	2-614	18
2-5130	7	1-5188TC	10	1-534	13	1-616	18
1-5130TC	7	1-5192	10	1-536	13	2-616	18
1-5134	7	2-5192	10	1-536TC	13	1-617	18
1-5134TC	7	1-5192TC	10	1-539	13	1-640	18
1-5138	7	1-5312	10	1-552	13	1-642	18
1-5138TC	7	1-5314	10	2-552	13	2-642	18
1-103-5145	7	1-5300	11	1-554	13	1-644	18
1-5144	8	2-5300	11	1-556	13	2-644	18
1-5144TC	8	1-5304	11	1-556TC	13	1-646	18
1-5146	8	1-5304	11	1-9105	14	2-646	18
1-5146TC	8	1-5304TC	11	1-9104	14	1-648	18
1-5160	8	1-5304LTC	11	2-9104	14	2-648	18
1-5162	8	1-5301	11	1-9100	14	1-650	18
1-271001	8	1-5302	11	2-9100	14	2-650	18
1-271000	8	2-5302	11	1-9102	14	1-106-0545	18
1-5220	8	1-5306	11	1-9103	14	1-660	18
1-5220TC	8	2-5306	11	2-9103	14	1-662	18
1-5222	8	1-5306TC	11	1-9101	14	2-662	18
1-5222TC	8	1-5306LTC	11	2-9101	14	1-664	18
2-181652	9	1-181462	11	1-996	15	2-664	18
1-5177	9	1-181466	11	1-9124	15	1-666	18
1-5179	9	2-181466	11	2-9124	15	1-6124	19
2-5179	9	1-512	12	1-9120	15	2-6124	19
		2-512	12				

INDEX

1-6124XL	19	1-7252	23	1-7166	25	1-25812	28
1-6126	19	1-7254	23	1-7167	25	1-24804	28
1-6125	19	1-714	23	2-7167	25	1-497230	29
1-24573	19	2-714	23	1-792	26	1-497231	29
1-24570	19	1-716	23	1-794	26	1-497232	29
1-24571	19	2-716	23	2-794	26	1-1699	29
1-24550	19	1-736	24	1-7118	26	1-24703	29
1-24572	19	2-736	24	1-7120	26	1-842	30
1-6118	20	1-738	24	1-7122	26	1-842TC	30
2-6118	20	2-738	24	1-7124	26	1-844	30
1-6118D	20	1-742	24	1-7126	26	2-844	30
1-172510	20	2-742	24	1-7128	26	1-844TC	30
1-6120	20	1-746	24	1-7130	26	1-844L	30
2-6120	20	2-746	24	1-7132	26	1-844LTC	30
1-6120L	20	1-744	24	1-7136	26	1-846	30
1-6120XL	20	2-744	24	2-7136	26	2-846	30
1-6123XL	20	1-748	24	1-7138	26	1-846TC	30
1-172500	20	2-748	24	2-7138	26	1-848	30
#2 - 17302	21	1-7170	24	1-7140	26	2-848	30
#3 - 17303	21	2-7170	24	1-7142	26	1-848TC	30
#4 - 17304	21	1-7180	24	2-7142	26	1-8130TC	30
#5 - 17305	21	2-7180	24	1-7144	26	1-814	30
#6 - 17306	21	1-1163710	24	1-7146	26	1-814TC	30
#7 - 17307	21	1-1163720	24	1-7148	26	2-814TC	30
1-6154	21	1-7172	24	1-16109	27	1-814ATC	30
1-6104	21	2-7172	24	1-16160	27	1-814LTC	30
1-6106	21	1-7182	24	2-16160	27	1-815	30
1-6107	21	2-7182	24	1-16110	27	2-815	30
1-19371	22	1-1163725	24	1-16162	27	1-815TC	30
1-6142	22	1-1163730	24	1-243261	27	2-815TC	30
2-6142	22	1-7260	25	1-163	27	1-815LTC	30
1-6144	22	1-7264	25	1-164	27	1-816	30
1-6145	22	1-7262	25	1-166	27	2-816	30
1-6312	22	1-7266	25	2-166	27	1-816TC	30
1-6314	22	1-770	25	1-168	27	2-816TC	30
1-724	23	1-7150	25	1-1610	27	1-816LTC	30
2-724	23	2-7150	25	1-1611	27	1-817	30
1-724D	23	1-7152	25	1-1612	27	2-817	30
1-726	23	2-7152	25	1-1628	27	1-817TC	30
2-726	23	1-7154	25	1-1644	28	2-817TC	30
1-726D	23	2-7154	25	1-1646	28	1-828	31
1-72	23	1-7156	25	1-1648	28	1-828TC	31
2-72	23	1-7157	25	1-7210	28	1-830	31
1-78	23	2-7157	25	1-7200	28	1-830TC	31
1-7250	23	1-7160	25	2-7200	28	1-833TC	31
1-74	23	2-7160	25	1-7202	28	1-834TC	31
2-74	23	1-7162	25	1-7206	28	1-835TC	31
1-710	23	2-7162	25	1-7212	28	1-850	31
1-74M	23	1-7164	25	1-25810	28	1-850TC	31
		2-7164	25			1-8750TCLH	31

INDEX

1-851	31	1-411	36	1-1814	41	1-19575	46
1-851TC	31	2-411	36	1-1835	41	1-19576	46
1-852TC	31	1-417	36	2-1835	41	1-19577	46
1-132-313TC	31	1-7504	37	1-1836	41	1-10310	46
1-822	31	2-7504	37	2-1836	41	2-10310	46
2-822	31	1-7506	37	1-181620	42	1-10312	46
1-896	32	2-7506	37	1-181622	42	1-11600	47
1-896TC	32	1-7501	37	1-181624	42	1-11608	47
1-87ATC	32	1-7500	37	1-394441	42	1-11614	47
1-86	32	2-7500	37	1-39442	42	1-11617	47
1-86TC	32	1-7514	37	1-181576	42	1-11602	47
1-84	32	1-7515	37	1-181510	42	1-11610	47
2-84	32	2-7515	37	1-181516	42	1-11616	47
1-84TC	32	1-28131	37	1-18450	43	1-11618	47
1-88	32	2-28131	37	1-18452	43	1-11364	47
1-88TC	32	1-28131N	37	1-18454	43	2-11364	47
2-8688TC	33	1-28132WB	37	1-181820	43	1-11362	47
1-870	33	1-28133	37	1-181828	43	1-11365	47
1-870TC	33	1-1718	37	2-181828	43	2-11365	47
1-181820	33	1-97590	38	1-181828TC	43	1-11366	47
1-181828	33	1-7604	38	1-181824	43	2-11366	47
2-181828	33	1-7600	38	1-181832	43	1-11367	47
1-181828TC	33	2-7600	38	1-181832TC	43	1-11368	47
1-181824	33	1-7640	38	1-181834	43	1-1112	48
1-181832	33	1-7614	38	1-181836	43	2-1112	48
1-181832TC	33	1-7610	38	1-181838	43	1-1113	48
1-181834	34	2-7610	38	2-181838	43	1-116	48
1-181836	34	1-18650	39	1-181839	43	1-112	48
1-181838	34	1-181202	39	1-18110	43	1-1114	48
2-181838	34	2-181202	39	1-18112	43	1-1116	48
1-181839	34	1-18833	39	1-18114	43	1-11430	49
1-8101TC	34	1-18866	39	1-18116	43	1-11420	49
1-8103TC	34	1-18867	40	1-5150027	44	2-11420	49
1-8104TC	34	2-18867	40	1-1872	44	1-11421	49
1-8105TC	34	1-18869	40	1-1874	44	1-1174	49
1-8106TC	34	1-18868	40	1-18499	44	1-1176	49
1-8107TC	34	1-18785	40	1-18500	44	2-1176	49
1-866TC	34	1-18786	40	1-18502	44	1-11114	49
1-869TC	34	1-18788	40	1-18504	44	1-11116	49
1-47	35	1-18780	40	1-18506	44	1-11122	49
2-47	35	1-18782	40	1-26782	45	2-11122	49
1-48	35	1-18784	40	1-19580	45	1-11180	50
2-48	35	1-357647	40	1-19582	45	1-11182	50
1-49	35	1-357645	40	1-19584	45	1-11184	50
2-49	35	1-357681TC	40	1-2104	45	1-1173	50
1-416	35	1-18827	41	1-2110	45	2-1173	50
2-416	35	1-18950	41	1-19570	46	1-11502	50
1-410	36	1-18951	41	1-19572	46	1-11530	50
2-410	36	1-1812	41	1-19573	46	1-11532	50
				1-19574	46	1-11534	50

INDEX

1-11536	50	1-275	57	1-19718	60	1-192100	67
1-11503	50	1-274	57	1-19720	60	1-19405	67
1-11538	50	1-276	57	1-157-6313	60	2-19405	67
1-11540	50	1-272	57	1-157-6316	60	1-19410	67
1-11542	50	1-2721	57	1-157-6319	60	2-19410	67
1-11544	50	1-2723	57	1-157-6323	60	1-19415	67
1-27213	51	1-2720	57	1-157-6325	60	1-19260	67
1-27214	51	1-2722	57	1-27700	61	1-19264	67
1-27215	51	1-49800	58	1-27702	61	1-19268	67
1-27216	51	1-27158TC	58	1-21638	61	2-19268	67
1-27217	51	1-27418	58	1-27703	61	1-19280	67
1-27218	51	1-27420	58	1-19804	61	1-19282	67
1-11250	51	1-27422	58	1-26190	62	1-19286	67
1-11252	51	1-27424	58	1-26192	62	1-19270	67
1-11254	51	1-27426	58	1-19808	62	1-19274	67
1-11256	51	1-27427	58	1-191230	62	1-19276	67
1-11260	51	1-27428	58	1-191232	62	1-19290	67
1-11112	52	1-27430	58	1-174800	63	1-19292	67
2-11112	52	1-27432	58	1-174801	63	2-19292	67
1-11700	52	1-27438	59	1-271390	63	1-19293	67
1-11200	52	1-27440	59	1-271391	63	1-19296	67
1-11202	52	1-27442	59	1-271392	63	2-19296	67
1-11204	52	1-27444	59	1-271394	63	1-19314	68
1-11206	52	1-27446	59	1-25454	63	2-19314	68
1-11208	52	1-27447	59	1-25456	63	1-19316	68
1-11210	52	1-27448	59	1-25458	63	2-19316	68
1-11212	52	1-27450	59	1-21626	64	1-19318	68
1-11216	52	1-27452	59	1-21626	64	2-19318	68
1-11102	53	1-27488	59	1-25386	64	1-181968	68
1-18102	53	1-27490	59	1-21627	64	1-191330	68
1-18100	53	1-27492	59	1-25385	64	1-191334	68
1-11138	54	1-27494	59	1-19844	64	1-2028	68
1-1106535C	54	1-27499	59	1-19724	64	1-2030	68
1-21153	54	1-27500	59	1-271295	64	1-3400	69
1-21154	54	1-27502	59	1-301000	65	1-3404	69
1-21156	54	1-27504	59	1-7620	65	1-3408	69
1-21158	54	1-27820	59	1-301420	65	1-3400D-2	69
1-21160	54	1-27821	59	1-301425	65	1-3404-2	69
1-2185	55	1-27822	59	1-28294	65	1-3408-2	69
1-2186	55	1-27824	59	1-28296	65	1-1072	69
1-2187	55	1-27826	59	1-28298	65	2-1072	69
1-18416	55	1-261410	60	1-3005	66	1-1074	69
1-261084	55	1-261412	60	1-3010	66	2-1074	69
1-19790	56	1-19704	60	1-3015	66	450205	70
1-19792	56	1-19706	60	1-3020	66	381-201	70
1-19793	56	1-19708	60	1-3040	66	381-202	70
1-11312	56	1-19710	60	1-3050	66	381-228	70
1-69846	56	1-19712	60	1-3055	66	381-269	70
1-27160	57	1-19714	60	1-3060	66	381-285	70
1-277	57	1-19716	60	1-302350	66	381-301	70

INDEX

T-01-01	70	2-6124	76	2-7504	77	SOP-2	78
T-01-02	70	1-644	76	1-5124	77	1-18867	78
T-01-03	70	2-644	76	2-5124	77	2-18867	78
T-01-04	70	1-7504	76	1-72	77	1-181202	78
T-01-05	70	2-7504	76	2-72	77	2-181202	78
T-01-07	70	1-5124	76	1-74	77	1-1835	78
T-01-09	70	2-5124	76	2-74	77	2-1835	78
465160	71	1-5182	76	1-736	77	1-181466	78
057150	71	2-5182	76	2-736	77	2-181466	78
KM47870	71	1-74	76	1-738	77	1-181828	78
181816	72	2-74	76	2-738	77	2-181828	78
510440	72	1-738	76	1-166	77	1-1310	78
510445	72	2-738	76	2-166	77	1-176400	78
KM71111	72	1-7172	76	1-47	77	CP634	78
KM71112	72	2-7172	76	2-47	77	AOP-1	78
KM71113	72	1-166	76	1-28131	77	AOP-2	78
KM35202	73	2-166	76	2-28131	77	1-27213	78
KM35210	73	1-47	76	CSP-1	77	1-27214	78
KM35212	73	2-47	76	CSP-2	77	1-27216	78
KM35214	73	1-28131	76	1-815	77	1-714	78
KM35220	73	2-28131	76	2-815	77	2-714	78
KM35250	73	SP-1	76	1-6124	77	1-1174	78
KM35280	74	SP-2	76	2-6124	77	1-11364	78
KM35292	74	1-815TC	76	1-7504	77	2-11364	78
KM35316	74	2-815TC	76	2-7504	77	1-11365	78
KM35286	74	1-6124	76	1-5124	77	2-11365	78
KM35298	74	2-6124	76	2-5124	77	KNEE-1	78
KM35322	74	1-7506	76	1-72	77	KNEE-2	78
KM35340	74	2-7506	76	2-72	77	1-19708	79
KM35364	74	1-5124	76	1-74	77	1-815	79
KM35380	74	2-5124	76	2-74	77	2-815	79
KM35392	74	1-5182	76	1-736	77	1-6124	79
KM35370	74	2-5182	76	2-736	77	2-6124	79
KM35382	74	1-74	76	1-738	77	1-5182	79
KM35394	74	2-74	76	2-738	77	2-5182	79
KM35396	74	1-736	76	1-47	77	1-74	79
KM35406	74	2-736	76	2-47	77	2-74	79
KM35418	74	1-7172	76	1-28131	77	1-22480	79
KM35430	74	2-7172	76	2-28131	77	1-11600	79
KM35450	74	1-47	76	FSP-1	77	1-11608	79
KM35454	74	2-47	76	FSP-2	77	1-47	79
KM35452	74	1-28131	76	1-84	78	2-47	79
KM35460	74	2-28131	76	2-84	78	1-1174	79
KM35484	75	GSP-1	76	1-18867	78	EAR-1	79
KM35490	75	GSP-2	76	2-18867	78	EAR-2	79
KM35494	75	1-816	77	1-1835	78	1-816	79
KM35496	75	2-816	77	2-1835	78	2-816	79
1-815	76	1-6124	77	1-181466	78	1-610	79
2-815	76	2-6124	77	2-181466	78	2-610	79
1-6124	76	1-7504	77	SOP-1	78	1-6124	79

INDEX

2-6124	79	2-11366	80	80.FR-12-04	91	005200	97
1-6118	79	1-11374	80	80.FR-12-05	91	001200	97
2-6118	79	1-410	80	80.FR-12-07	91	001204	97
1-7515	79	2-410	80	80.B-16-030	91	002200	97
2-7515	79	ETP-1	80	80.B-16-040	91	002204	97
1-5124	79	ETP-2	80	80.B-16-050	91	003200	97
2-5124	79	1-171026TC	81	80.B-20-050	91	003204	97
1-5181	79	1-394441	81	80.B-16-120	91	004200	97
2-5181	79	357660	81	80.B-16-100	91	004204	97
1-72	79	357661	81	80.D-10-159	91	BFWP300	98
2-72	79	1-357647	81	80.D-16-103	91	BFW5500	98
1-7180	79	1-47	81	80.TB-2100	91	BFW9030	98
2-7180	79	2-47	81	80.TN-2121	91	BFW9030B	98
1-1173	79	1-5312	81	80.KSLSR-8501	92	BFW9030M	98
2-1173	79	1-6118	81	80.KSLSR-8512	92	DE571	99
1-47	79	2-6118	81	80.3FRIN556120	92	ES140	99
2-47	79	1-1112	81	80.3FRIN556125	92	ES108	99
ESP-1	79	2-1112	81	80.T-SAMP-RMV	92	ES140	99
ESP-2	79	CP1038S	81	80.T-08RMV	92	DE571	99
1-817	80	MIP-1	81	80.T-32RMV	92	ES108	99
2-817	80	MIP-2	81	80.11000-000	93	78000-HP	101
1-8102TC	80	CP614-DK-1	83	80.11000-002	93	78001-HP	101
1-616	80	CP614-DK-2	83	80.11000-003	93	78002-HP	101
2-616	80	CP614-DK-3	83	80.11000-500	93	73000	101
1-5190	80	CP614	83	80.11000-506	93	4005B	101
1-1612	80	CP778-DK-1	84	80.11000-507	93	4005S	101
1-7144	80	CP778-DK-2	84	80.11000-112	93	4012B	101
1-410	80	CP778-DK-3	84	80.11000-400	93	4012B-2	101
2-410	80	CP778	84	80.3FSZB790111	93	4005S-2	101
1-411	80	CP634	85	80.K1492-20	94	4005B-2	101
2-411	80	CP634-DK-5	85	80.K4387-20	94	06610	102
ELP-1	80	CP1038D2	85	80.K1811-20	94	1018/Rev2	103
ELP-2	80	CP1038D1	86	528135	95	1024/Rev2	103
1-817	80	CP1038D1-P-1S/L	86	528235	95	1030/Rev 2	103
2-817	80	CP1038D1-P-3S/L	86	525980	95	1036/Rev 2	103
1-616	80	CP1038S	86	537061	96	1048/Rev 2*	103
2-616	80	CP471D1-P-1S/L	87	537081	96	1060/Rev 2*	103
1-5190	80	CP471D1-P-3S/L	87	137061	96	1072/Rev 2**	103
1-5134	80	CP471D1	87	137081	96		
1-5222	80	CP471S-DK-2	88	237061	96		
1-748	80	CP471S-DK-3	88	237081	96		
2-748	80	CP471S-DK-4	88	237117	96		
1-1612	80	CP471S-DK-5	88	337081	96		
1-7126	80	CP471S	88	337085	96		
1-7144	80	REC128	89	337117	96		
1-1648	80	ACG15	90	437081	96		
1-7600	80	ACG17	90	437085	96		
2-7600	80	80.M-08-005	91	437117	96		
1-27920	80	80.M-12-020	91				
1-11366	80	80.FR-12-03	91				



RICA Animal Health is a Division of RICA Surgical

**9207 Ivanhoe Street
Schiller Park, IL 60176**

**847-671-5600
1-800-889-3218
Fax: 847-671-5690**

www.RICAanimalhealth.com



RICA Animal Health is a Pacific Surgical Company

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.RICAanimalhealth.com

* Universal Wing Adapter Assembly is Patent Pending.

©2025 RICA Animal Health
PRINTED IN THE U.S.A.

FDA Registered Facility

