

# SAFE & EFFECTIVE WARMTH FOR PETS

*Prevent Hypothermia and Keep Patients Warm with RICA's Patented\* Thermal Controlled Products*

- ⇒ **Safe and consistent temperature of 98° - 103° F.**
- ⇒ **More durable.** Built with heavy duty ABS medical grade plastic.
- ⇒ **Portable and user-friendly.** Simple and quick to move and use anywhere in the veterinary hospital. Plug into any electrical outlet.
- ⇒ **Supported by RICA's 3-year conditional warranty.** *See website for further details.*

## Surgery/Recovery/Dental Thermal Pads



- ▶ Available in four useful and convenient sizes.  
12" x 18"                      20" x 20"  
20" x 28"                      15" x 45"
- ▶ Fits in all size cages. Each pad comes with a 6-foot cord with added strain relief protection to prevent the cord from being pulled out.
- ▶ Excellent to use for birds, small mammals, all dental procedures, post-op recovery, etc.
- ▶ Easy to clean and sanitize. **Water-resistant.**

## V-Top Surgical Table Thermal Pad



- ▶ Use on a V-Top or flat surgery table.
- ▶ Consistent temperature and safe heat distribution throughout the entire pad.
- ▶ Versatile. Able to easily move and use elsewhere in the clinic.
- ▶ Measures 18.5" wide x 41" long.

**Contact your RICA Surgical Products Veterinary Distributor to learn more about these innovative and cost-effective Thermal Controlled Products.**

\*U.S. Patent # 5,371,340



RICA® Surgical Products, Inc.  
1-800-889-3218  
www.thermalpads4pets.com