




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*


Item CP539060 Universal Wing Adapter Assembly*




The Universal Wing Adapter Assembly includes one Universal Wing Adapter, one spring, two nuts, two 5.9 mm screws, and two 7.8 mm screws (for use in thicker aluminum and plastic sterilization trays). This item has a limited lifetime and will need to be replaced occasionally due to wear and tear.

*Universal Wing Adapter Assembly is Patent Pending


Item CP530000 Family of Silicone Holders



The **CP530000 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP530001 1.5", One 3mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP530002 2", Two 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP530003 3", Three 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530004 4", Three 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530005 5", Four 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530006 6", Five 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530007 6", Six 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530008 7", Six 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP530200 Family of Silicone Holders



The **CP530200 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.

©2020

FORM # SBUSIPS-530200-0420

CP530200 Family



Item CP530201 1.5", One 5mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP530202 1.5", Two 5mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP530203 2", Three 5mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP530204 3", Four 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530205 3", Five 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP530206 4", Seven 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP530207 5", Eight 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP530208 5", Nine 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530209 6", Ten 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530210 6", Eleven 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP530211 7", Twelve 5mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP537000 Family of Silicone Holders



The **CP537000 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP537001 1.5", One 3mm-8mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537002 1.5", Two 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537003 2", Two 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly

CP537000 Family



Item CP537004 2", Three 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537005 3", Three 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537006 3", Four 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537007 4", Five 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537008 4", Six 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537009 5", Six 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537010 5", Seven 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537011 6", Eight 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537012 6", Nine 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537013 7", Ten 3mm-8mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.




Universal Scope & Instrument Protection System

from RICA Surgical Products


SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP537100 Family of Silicone Holders



The **CP537100 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.

CP537100 Family



Item CP537101 1.5", One 8mm-15mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537102 2", Two 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537103 3", Three 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537104 4", Four 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537105 4", Five 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537106 5", Five 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537107 5", Six 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537108 6", Seven 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537109 7", Eight 8mm-15mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.

©2020

FORM # SBUSIPS-537100-0420




Universal Scope & Instrument Protection System


from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP537200 Family of Silicone Holders



The **CP537200 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP537201 1.5", One 5mm & 10mm Double Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537202 2", Two 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP537203 3", Three 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537204 3", Four 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537205 4", Four 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537206 4", Five 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537207 5", Six 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537208 5", Seven 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537209 6", Eight 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537210 7", Nine 5mm & 10mm Double Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP537300 Family of Silicone Holders



The **CP537300 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP537301 1.5", One 3mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP537302 1.5", Two 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537303 2", Two 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly


CP537300 Family



Item CP537304 2", Three 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537305 3", Four 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537306 3", Five 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537307 4", Five 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537308 4", Six 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537309 5", Eight 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537310 5", Nine 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537311 6", Ten 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537312 6", Eleven 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537313 7", Twelve 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Item CP537400 Family of Silicone Holders




The **CP537400 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP537401 1.5", One 10mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537402 2", Two 10mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP537403 3", Three 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537404 4", Four 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537405 5", Five 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537406 5", Six 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537407 6", Six 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537408 6", Seven 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537409 7", Eight 10mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.




Universal Scope & Instrument Protection System

from RICA Surgical Products

SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*


Item CP537500 Family of Silicone Holders



The **CP537500 Family** of silicone holders provides multiple configurations to secure and protect your scopes and instruments, minimize scope and instrument movement, and maximize sterilant flow and contact.



Item CP537501 1.5", One 3mm Opening, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP537502 1.5", Two 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP537503 2", Three 3mm Openings, Silicone Instrument Bar with Radel Bar and Universal Wing Adapter Assembly

CP537500 Family



Item CP537504 3", Four 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537505 4", Five 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537506 4", Six 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP537507 5", Seven 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537508 5", Eight 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537509 6", Nine 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537510 6", Ten 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP537511 7", Twelve 3mm Openings, Silicone Instrument Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.



Universal Scope & Instrument Protection System

from RICA Surgical Products


SteriBest® Universal Scope & Instrument Protection System has been designed to hold a wide variety of scopes and surgical instruments in virtually all trays and baskets used for sterilization and transport utilizing the Universal Wing Adapter Assembly.*

Silicone Scope Holders

The **SteriBest® Family of Silicone Scope Holders** provide multiple configurations to secure and protect your scopes, minimize scope movement, and maximize sterilant flow and contact.




Item CP535000 1.5", One 15mm Scope Holder, Open End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP535001 2", One 15mm Scope Holder, Open End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP535100 1.5", One 8mm Scope Holder, Closed End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP535101 2", One 8mm Scope Holder, Closed End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP535200 1.5", One 12mm Scope Holder, Closed End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP535201 2", One 12mm Scope Holder, Closed End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP535300 1.5", One 12mm Scope Holder, Open End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly




Item CP535301 2", One 12mm Scope Holder, Open End, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly



Scope Holders



Item CP535510 3", One 15mm Scope Holder, Open End, Silicone Bar and One 8mm Scope Holder, Closed End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP535511 3", One 12mm Scope Holder, Closed End, Silicone Bar and One 15mm Scope Holder, Open End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP535512 3", One 12mm Scope Holder, Open End, Silicone Bar and One 8mm Scope Holder, Closed End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP535513 3", One 12mm Scope Holder, Closed End, Silicone Bar and One 12mm Scope Holder, Open End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP535520 4", One 8mm Scope Holder, Closed End, Silicone Bar and One 15mm Scope Holder, Open End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies




Item CP535521 4", One 15mm Scope Holder, Open End, Silicone Bar and One 12mm Scope Holder, Closed End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP535522 4", One 8mm Scope Holder, Closed End, Silicone Bar and One 12mm Scope Holder, Open End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP535523 4", One 12mm Scope Holder, Open End, Silicone Bar and One 12mm Scope Holder, Closed End, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies



Item CP251300 2", Three 5.5mm Openings, Silicone Bar with Radel Bar and Universal Wing Adapter Assembly



Item CP461500 4", Five 5.5mm Openings, Silicone Bar with Radel Bar and 2 Universal Wing Adapter Assemblies

*Universal Wing Adapter Assembly is Patent Pending.

SteriBest® is a Registered Trademark of RICA Surgical Products, Inc., an FDA Registered Facility since 1984. Registration No. 1419732
RICA Surgical Products, Inc. Schiller Park, IL 60176 (800) 889-3218 www.ricasurgical.com
PRINTED IN THE U.S.A.