



RICA Surgical Products

# Weck® Horizon® Ligation Clip System

## Precise Easy-Loading Titanium Clips

Weck's Horizon System provides veterinary surgeons all of the benefits and time-savings associated with the Hemoclip System, with two additional benefits:

- Easier to load clip, since the design of the Horizon Cartridge creates a near frictionless one-hand load.
- An additional clip size...Small/Wide...for more precise ligation of the blood vessel with the appropriate size clip.



### The Horizon Clip

The pre-formed chevron shape of the Horizon clip engulfs tissue with precise tip-to-tip closure. The heart-shaped wire is designed to give each clip a firm grip on vessels, while the triangulated cross-section of the clip leg maximizes surface-to-surface contact between clip and jaw, practically eliminating clip fallout.

### The Horizon Cartridge

The design of the Horizon Cartridge creates a near frictionless one-hand load simplifying the loading of the clip into the applicator. Horizon Cartridges are color-coded to match applicators and the adhesive backing allows convenient placement of the cartridge to any sterile area.

*The Horizon System is the premium Weck titanium clip ligation platform, providing veterinary surgeons a strong, secure clip coupled with an efficient, convenient easy-load cartridge system.*



See reverse side for current pricing and packaging of the sterile, titanium Horizon clips.

# Weck® Horizon® Ligation Clip System 2022 Veterinary Price List

| LIST # | HORIZON CLIP SIZE | CLIPS/<br>CARTRIDGE | CARTRIDGES/<br>BOX | CLIPS/<br>BOX | VETERINARY PRICE<br>PER BOX | PER CLIP |
|--------|-------------------|---------------------|--------------------|---------------|-----------------------------|----------|
|--------|-------------------|---------------------|--------------------|---------------|-----------------------------|----------|

## Horizon Ligation Clip System – Titanium Clips


|        |                      |   |    |     |          |        |
|--------|----------------------|---|----|-----|----------|--------|
| 005200 | MicroClip™ (white)   | 6 | 30 | 180 | \$213.50 | \$1.19 |
| 001200 | Small (yellow)       | 6 | 30 | 180 | \$184.25 | \$1.02 |
| 001201 | Small-Wide (red)     | 6 | 30 | 180 | \$201.95 | \$1.12 |
| 002200 | Medium (blue)        | 6 | 30 | 180 | \$217.80 | \$1.21 |
| 003200 | Medium/Large (green) | 6 | 20 | 120 | \$186.15 | \$1.55 |
| 004200 | Large (orange)       | 6 | 20 | 120 | \$249.60 | \$2.08 |

|        |                      |    |    |     |          |        |
|--------|----------------------|----|----|-----|----------|--------|
| 001204 | Small (yellow)       | 24 | 25 | 600 | \$560.75 | \$0.93 |
| 001205 | Small-Wide (red)     | 24 | 25 | 600 | \$612.30 | \$1.02 |
| 002204 | Medium (blue)        | 24 | 25 | 600 | \$661.85 | \$1.10 |
| 003204 | Medium/Large (green) | 24 | 15 | 360 | \$507.50 | \$1.41 |
| 004204 | Large (orange)       | 24 | 15 | 360 | \$680.60 | \$1.89 |

| LIST # | HORIZON APPLIERS     | LENGTH | TYPE               | PRICE    |
|--------|----------------------|--------|--------------------|----------|
| 537061 | MicroClip Applier    | 6.0"   | Curved (20°)       | \$172.50 |
| 537081 | MicroClip Applier    | 8.0"   | Curved             | \$172.50 |
| 137061 | Small Applier        | 6.0"   | Curved             | \$172.50 |
| 137081 | Small Applier        | 8.0"   | Curved             | \$172.50 |
| 137062 | Small/Wide Applier   | 6.0"   | Curved             | \$172.50 |
| 137082 | Small/Wide Applier   | 8.0"   | Curved             | \$172.50 |
| 237061 | Medium Applier       | 6.0"   | Curved             | \$172.50 |
| 237081 | Medium Applier       | 8.0"   | Curved             | \$172.50 |
| 237117 | Medium Applier       | 11.0"  | Right Angled (70°) | \$172.50 |
| 337081 | Medium/Large Applier | 8.0"   | Curved             | \$172.50 |
| 337085 | Medium/Large Applier | 8.0"   | Angled (55°)       | \$172.50 |
| 337117 | Medium/Large Applier | 11.0"  | Right Angled       | \$172.50 |
| 437081 | Large Applier        | 8.0"   | Curved             | \$172.50 |
| 437085 | Large Applier        | 8.0"   | Angled             | \$172.50 |
| 437117 | Large Applier        | 11.0"  | Right Angled       | \$172.50 |

## DIMENSIONS OF WECK HORIZON CLIPS

| HORIZON DESCRIPTION       | WIDTH  |             | HEIGHT |             | CLOSED CLIP LENGTH |             |
|---------------------------|--------|-------------|--------|-------------|--------------------|-------------|
|                           | INCHES | MILLIMETERS | INCHES | MILLIMETERS | INCHES             | MILLIMETERS |
| MicroClip Horizon (white) | 0.075  | 1.91        | 0.075  | 1.91        | 0.088              | 2.24        |
| Small Horizon (yellow)    | 0.078  | 1.98        | 0.113  | 2.87        | 0.143              | 3.63        |
| Small-Wide Horizon (red)  | 0.082  | 2.08        | 0.113  | 2.87        | 0.145              | 3.68        |
| Medium Horizon (blue)     | 0.119  | 3.02        | 0.184  | 4.67        | 0.232              | 5.89        |
| Med-Large Horizon (green) | 0.226  | 5.74        | 0.284  | 7.21        | 0.353              | 8.97        |
| Large Horizon (orange)    | 0.307  | 7.80        | 0.396  | 10.07       | 0.488              | 12.38       |




NOTE: Width and height are taken from **INSIDE** the clip. For informational use only, all dimensions are **APPROXIMATE**.

**To order, contact your authorized RICA Veterinary Distributor for current pricing, call RICA's office in Schiller Park, Illinois at 1-800-889-3218 or visit [www.ricasurgical.com](http://www.ricasurgical.com).**



**RICA Surgical Products**

**RICA Surgical Products, Inc.**  
9207 Ivanhoe Street  
Schiller Park, IL 60176  
(800) 889-3218



PRINTED IN USA  
FORM # RSP-WHLCS-0122